

Datasheet: HCA108

Description:	RECOMBINANT HUMAN IgG2 LAMBDA
Name:	HUMAN IgG2 LAMBDA
Format:	Purified
Product Type:	Recombinant Protein
Clone:	AbD264_hIgG2
Isotype:	IgG2 Lambda
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			
Immunoprecipitation			▪	
Western Blotting			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Human IgG2 antibody selected from the HuCAL® phage display library and expressed in a human cell line, HKB-11. This antibody is supplied liquid.
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Thiomersal
Approx. Protein Concentrations	Total protein concentration 0.5 mg/ml
Immunogen	Green fluorescent protein
External Database Links	UniProt: C6KXN3 Related reagents P01859 Related reagents

Entrez Gene:[3501](#) IGHG2 [Related reagents](#)

Product Information	<p>Recombinant Human IgG2 lambda contains a lambda light chain. The heavy chain belongs to the G2m(..), previously designated G2m(n-), allotype.</p> <p>Recombinant Human IgG2 lambda is a recombinant antibody with specificity for Green Fluorescent Protein. It has no known reactivity with mammalian proteins or other antigens. It is therefore recommended as a control reagent for assays where this immunoglobulin is measured, or when using other human antibodies of the same isotype and subclass. It is recommended that this reagent is used at the same concentration as the test reagent.</p>
References	<p>1. Madan, B. <i>et al.</i> (2016) USP6 oncogene promotes Wnt signaling by deubiquitylating Frizzleds. Proc Natl Acad Sci U S A. May 9. pii: 201605691. [Epub ahead of print]</p>
Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Shelf Life	<p>12 months from date of despatch.</p>
Acknowledgements	<p>Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.</p>
Health And Safety Information	<p>Material Safety Datasheet Documentation #10094 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf</p>
Licensed Use	<p>For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.</p>
Regulatory	<p>For research purposes only</p>
Technical Advice	<p>Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual</p>

Related Products

Recommended Secondary Antibodies

Mouse Anti Human IgG (Fc) (CH2 DOMAIN) (MCA647...)	HRP
Goat Anti Human IgG F(ab') ₂ (STAR126...)	Purified , Alk. Phos. , Biotin , FITC , HRP , TRITC
Goat Anti Human IgG F(ab') ₂ (0500-0099...)	Alk. Phos. , HRP

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

"M304967:170420"

Printed on 17 May 2017

