

Datasheet: HCA144

Description:	HUMAN ANTI HUMAN CD127
Specificity:	CD127
Other names:	IL-7R
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD11590
Isotype:	HuCAL Fab bivalent
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	▪			1/10 - 1/200
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting		▪		
Functional Assays			▪	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)HCA144 has been demonstrated to function in Flow Cytometry, however Bio-Rad recommend the use of clone 11592 (HCA145) which demonstrates stronger staining.

Target Species	Human
Product Form	A bivalent human recombinant Fab selected from the HuCAL® phage display library, expressed in <i>E. coli</i> . This Fab fragment is dimerized via a helix-turn-helix motif. The antibody is tagged with a DYKDDDDK tag and a double extended Strep-tag at the C-terminus of the antibody heavy chain. This antibody is supplied lyophilised.
Reconstitution	Reconstitute with 1.0ml distilled water
Preparation	StrepTactin affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide (NaN ₃)
Stabilisers	1% Bovine Serum Albumin

Approx. Protein Concentrations	Total protein concentration 0.1mg/ml
Immunogen	Recombinant Human IL-7 R alpha/CD127/Fc Chimera.
External Database Links	<p>UniProt: P16871 Related reagents</p> <p>Entrez Gene: 3575 IL7R Related reagents</p>
Specificity	Human anti Human CD127 antibody, clone 11590 recognizes the alpha chain of the human interleukin-7 receptor (IL-7R), which is also known as CD127. CD127 is a ~65-75 kDa transmembrane glycoprotein that associates with the interleukin-2 receptor gamma chain to form the functional high affinity IL-7 receptor. CD127 is expressed on lymphoid precursors in the bone marrow, pro-B cells, a subpopulation of thymocytes, T-cells and monocytes. CD127 is the specific receptor for IL-7 and the IL-7 receptor plays a critical role in lymphoid development.
ELISA	Due to the presence of bovine serum albumin (BSA), this antibody is unsuitable for use as a capture reagent in sandwich ELISA applications. This product is also available without BSA, please enquire .
Histology Positive Control Tissue	Human tonsil
Storage	Prior to reconstitution store at +4°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of reconstitution.
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10041 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
Licensed Use	For in vitro research purposes only, unless otherwise specified in writing by Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Secondary Antibodies

Rabbit Anti Synthetic Peptide DYKDDDDK TAG (AHP1074...) [HRP](#)

Mouse Anti Synthetic Peptide STREP-TAG IMMO (MCA2488...) [Purified](#)

Recommended Useful Reagents

[HISTAR DETECTION SYSTEM \(STAR3000B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M304855:170420'

Printed on 13 Nov 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)