

Datasheet: HCA169G

Description:	HUMAN ANTI TRASTUZUMAB
Specificity:	TRASTUZUMAB
Other names:	HERCEPTIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AbD18141
Isotype:	HuCAL Fab monovalent
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form

A monovalent human recombinant Fab (lambda light chain) selected from the HuCAL® phage display library, expressed in E. coli. The antibody is tagged with a DYKDDDDK tag and a HIS-tag (HHHHHH) at the C-terminus of the antibody heavy chain. This antibody is supplied as a liquid.

Preparation

Metal chelate affinity chromatography

Buffer Solution

Phosphate buffered saline

Preservative Stabilisers

0.01% Thiomersal

Approx. Protein Concentrations

Antibody concentration 1.0 mg/ml

Immunogen

Trastuzumab

Specificity

Human Anti-Trastuzumab Antibody, clone AbD18141, is a paratope specific, high affinity, anti-idiotypic antibody that specifically recognizes the monoclonal antibody drug trastuzumab and inhibits it binding to its target. The antibody can be used to measure the levels of trastuzumab and biosimilar products in bioanalytical assays.

Clone AbD18141 is a fully human recombinant monoclonal antibody in monovalent Fab format and is suitable as a capture antibody paired with antibody clone AbD16712_hlgG1 ([HCA176P](#)) in a

pharmacokinetic (PK) bridging ELISA, to measure free drug in preclinical research and clinical trials using patient sera.

Trastuzumab, also known as Herceptin, is a drug used in the treatment of HER2 positive breast cancer and other HER2 over-expressing cancers including HER2-positive metastatic cancers of the gastrointestinal tract. Trastuzumab binds to the HER2 (or c-erbB2) proto-oncogene, an EGF receptor-like protein found on 20-30% of breast cancer cells. The binding leads to antibody mediated (complement mediated) killing of the HER2 positive cells.

[View a summary of all anti-trastuzumab antibodies.](#)

Affinity	The monovalent intrinsic affinity of this antibody was measured as $K_D=0.07$ nM by real time, label-free molecular interaction analysis on immobilized trastuzumab
ELISA	This product may be used as a capture antibody in a sandwich ELISA together with HCA176P as the detection reagent. Protocol: PK bridging ELISA to measure free drug
References	1. Bults, P. <i>et al.</i> (2016) LC-MS/MS-Based Monitoring of In Vivo Protein Biotransformation: Quantitative Determination of Trastuzumab and Its Deamidation Products in Human Plasma. Anal Chem. 88 (3): 1871-7.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Acknowledgements	Sold under license of U.S. Patents 6753136, 7785859 and 8273688 and corresponding patents. This antibody was developed by Bio-Rad, Zeppelinstr. 4, 82178 Puchheim, Germany. Herceptin is a registered trademark of Genentech, Inc.
Health And Safety Information	Material Safety Datasheet documentation #10094 available at: 10094: https://www.bio-rad-antibodies.com/uploads/MSDS/10094.pdf
Licensed Use	For in vitro research purposes and for commercial applications for the provision of in vitro testing services to support preclinical and clinical drug development. Any re-sale in any form or any other commercial application needs a written agreement with Bio-Rad.
Regulatory	For research purposes only
Technical Advice	Recommended protocols and further information about HuCAL recombinant antibody technology can be found in the HuCAL Antibodies Technical Manual

Related Products

Recommended Secondary Antibodies

Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA5995...) [HRP](#)

Rabbit Anti Synthetic Peptide DYKDDDDK TAG (AHP1074...) [HRP](#)

Goat Anti Synthetic Peptide HISTIDINE TAG (AHP1656...) [Biotin](#)
Rat Anti Synthetic Peptide DYKDDDDK TAG (MCA4764...) [HRP](#)
Mouse Anti Synthetic Peptide HISTIDINE TAG (MCA1396...) [Alk. Phos.](#), [Biotin](#)

Recommended Negative Controls

[HuCAL Fab-FH NEGATIVE CONTROL \(HCA045\)](#)

Recommended Useful Reagents

[HUMAN ANTI TRASTUZUMAB \(HCA176\)](#)
[HUMAN ANTI TRASTUZUMAB \(HCA168\)](#)
[HUMAN ANTI TRASTUZUMAB \(HCA177\)](#)
[HISPEC ASSAY DILUENT \(BUF049A\)](#)
[HUMAN ANTI TRASTUZUMAB \(HCA166\)](#)
[HUMAN ANTI TRASTUZUMAB \(HCA167\)](#)
[MOUSE ANTI HUMAN IgG \(Fc\) CH2 DOMAIN:HRP \(MCA647P\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

'M311625:171213'

Printed on 30 Apr 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)