

Datasheet: MCA1197B

Description:	MOUSE ANTI HUMAN CD7:Biotin
Specificity:	CD7
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	LT7
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative	0.09% Sodium Azide
Stabilisers	1% Bovine Serum Albumin
Approx. Protein Concentrations	IgG concentration 0.1mg/ml
Immunogen	Normal human peripheral blood lymphocytes.
External Database Links	UniProt: P09564 Related reagents Entrez Gene: 924 CD7 Related reagents
Fusion Partners	Spleen cells from immunized BALB/c mice were fused with cells of the X63.653 mouse myeloma

cell line.

Specificity	Mouse anti Human CD7 antibody, clone LT7 recognizes the human T-cell antigen CD7, also known as TP41 or CD7. CD7 is a single chain cell surface glycoprotein of ~40 kDa, expressed by T cells, NK cells and thymocytes. CD7 is also expressed by tumor cells in some cases of acute myeloid leukaemia.
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	1. Haynes, B.F. <i>et al.</i> (1989) Ontogeny of T-cell precursors: a model for the initial stages of human T-cell development. Immunol Today. 10 (3): 87-91. 2. Aruffo, A. & Seed, B. (1987) Molecular cloning of two CD7 (T-cell leukemia antigen) cDNAs by a COS cell expression system. EMBO J. 6 (11): 3313-6.
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from the date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M304496:170407'

Printed on 01 Aug 2018