

Datasheet: MCA1117

Description:	MOUSE ANTI HUMAN INTERFERON ALPHA
Specificity:	IFN ALPHA
Other names:	INTERFERON ALPHA
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	C10F5
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA	■			
Immunoprecipitation			■	
Western Blotting			■	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant Human Interferon alpha -A.
External Database Links	UniProt: P01562 Related reagents

[P01568](#) [Related reagents](#)

[P01563](#) [Related reagents](#)

Entrez Gene:

[3447](#) IFNA13 [Related reagents](#)

[3452](#) IFNA21 [Related reagents](#)

[3440](#) IFNA2 [Related reagents](#)

Fusion Partners Spleen cells from immunized BALB/c mice were fused with cells of the X63. Ag8.653 myeloma cell line.

Specificity **Mouse anti Human interferon alpha antibody, clone C10F5** recognizes an epitope common to IFN-alpha subtypes A, F, D, I1 and P. Strongest reaction in ELISA is with subtypes alpha F and alpha D. It has an affinity constant of 2.4×10^{-10} M.

Mouse anti Human interferon alpha antibody, clone C10F5 has been used successfully for the immunohistochemical detection of interferon- α producing cells on cryostat sections in a human suction blister model ([Reed et al. 2004](#)).

References

1. Reed, J.R. *et al.* (2004) Telomere erosion in memory T cells induced by telomerase inhibition at the site of antigenic challenge *in vivo*. [J Exp Med. 199 \(10\): 1433-43.](#)
2. van der Laan, A. M. (2013) Adverse left ventricular remodelling after acute myocardial infarction in patients is associated with attenuated type I interferon signalling In: [Monocytes in ischemic heart disease 's-Hertogenbosch: Boxpress](#)

Storage Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 18 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...) [RPE](#)

Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)

Rabbit Anti Mouse IgG (STAR9...) [FITC](#)

Goat Anti Mouse IgG (STAR77...) [HRP](#)

Rabbit Anti Mouse IgG (STAR12...) [RPE](#)

Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

Rabbit Anti Mouse IgG (STAR8...) [DyLight®800](#)

Goat Anti Mouse IgG (STAR70...) [FITC](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Human Anti Mouse IgG1 (HCA036...) [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight@488](#), [DyLight@549](#),
[DyLight@649](#), [DyLight@680](#), [DyLight@800](#),
[FITC](#), [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M319000:180719'

Printed on 01 Aug 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)