

Datasheet: MCA1269PE

Description:	MOUSE ANTI HUMAN CD100:RPE
Specificity:	CD100
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	A8
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Species Cross Reactivity	Reacts with: Cynomolgus monkey, Rhesus Monkey N.B. Antibody reactivity and working conditions may vary between species.		
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
Reconstitution	Reconstitute with 1 ml distilled water.		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant.		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1% Bovine Serum Albumin 5% Sucrose		
Immunogen	Human PHA activated lymphocytes.		
External Database Links	UniProt:		

[Q92854](#) [Related reagents](#)

Entrez Gene:

[10507](#) SEMA4D [Related reagents](#)

Synonyms	C9orf164, CD100, SEMAJ
-----------------	------------------------

Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS1 myeloma cell line.
------------------------	---------------------------------------------------------------------------------------------------

Specificity	Mouse anti Human CD100 antibody, clone A8 recognizes human Sepaphorin 4D, also known as CD100, A8, BB18 or GR3. CD100 is a 841 amino acid ~150 kDa homodimeric single pass type 1 transmembrane cell surface glycoprotein acting as a receptor for some members of the plexin family. CD100 is expressed at low levels by peripheral blood T cells, and at higher density on activated T cells and B cells, it is also expressed at high levels in fetal kidney with lower levels in fetal brain, liver and lung (Maestrini <i>et al.</i> 1996).
--------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
-----------------------	-------------------------------------------------------------------------------------

References	<ol style="list-style-type: none">1. Aversa, G.G. <i>et al.</i> (1989) Use of monoclonal antibodies to study <i>in vivo</i> and <i>in vitro</i>-activated lymphocytes. Transplant Proc. 21 (1 Pt 1): 349-50.2. Van Vliet, S. J. <i>et al.</i> (1995) Species cross reactivity (human-monkey-bovine) of the adhesion structure section mAb. Leucocyte Typing V. Oxford University Press p1607-8.3. Mizrahi, S. <i>et al.</i> (2007) CD100 on NK cells enhance IFNgamma secretion and killing of target cells expressing CD72. PLoS One. 2 (9): e818.4. Yoshino, N. <i>et al.</i> (2000) Upgrading of flow cytometric analysis for absolute counts, cytokines and other antigenic molecules of cynomolgus monkeys (<i>Macaca fascicularis</i>) by using anti-human cross-reactive antibodies. Exp Anim. 49 (2): 97-110.
-------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Storage	Prior to reconstitution store at +4°C. Following reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Shelf Life	12 months from date of reconstitution.
-------------------	----------------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10075 available at: 10075: https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf
--------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA928PE\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com
'M300580:170106'

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

Printed on 14 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)