

Datasheet: MCA1474T

Description:	RAT ANTI MOUSE THY-1.2
Specificity:	CD90
Other names:	THY1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	FF-10
Isotype:	IgG2c
Quantity:	25 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/50 - 1/100
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA			■	
Immunoprecipitation	■			
Western Blotting	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	EL-4 cells.

**External Database
Links**

UniProt:

[P01831](#) [Related reagents](#)

Entrez Gene:

[21838](#) Thy1 [Related reagents](#)

Synonyms

Thy-1

Fusion Partners

Spleen cells from immunised Sprague-Dawley rats were fused with cells from the mouse SP2/0 - Ag14 myeloma cell line.

Specificity

Rat anti Mouse THY-1.2 antibody, clone FF-10 recognizes the murine Thy1.2 cell surface antigen, and binds to T cell lines (EL-4, 3A9), T lymphocytes and the WEH1-3 cell line. 3T3 fibroblasts are also stained by this antibody.

Rat anti Mouse THY-1.2 antibody, clone FF-10 strongly induces homotypic adhesion of Thy1.2 positive cell lines.

Western Blotting

MCA1474T recognises a 30kD glycoprotein in immunoprecipitation and Western blotting assays.

References

1. Fan, X. *et al.* (1996) ICAM-1, -2, and LFA-1 independent homotypic T cell aggregation induced by a novel activating monoclonal antibody targeting the murine Thy-1 molecule. [Inflammation 20: 401.](#)
2. Venugopalan P *et al.* (2016) Transplanted neurons integrate into adult retinas and respond to light. [Nat Commun. 7: 10472.](#)
3. Szczerbal, I. *et al.* (2009) The spatial repositioning of adipogenesis genes is correlated with their expression status in a porcine mesenchymal stem cell adipogenesis model system. [Chromosoma. 118 \(5\): 647-63.](#)

Storage

Store at +4°C for one month or at -20°C for longer.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Rat IgG (STAR16...)	DyLight®800
Goat Anti Rat IgG (STAR73...)	RPE
Rabbit Anti Rat IgG (STAR21...)	HRP
Rabbit Anti Rat IgG (STAR17...)	FITC
Goat Anti Rat IgG (MOUSE ADSORBED) (STAR71...)	DyLight®549 , DyLight®649 , DyLight®800

Goat Anti Rat IgG (STAR131...)

[Alk. Phos.](#), [Biotin](#)

Goat Anti Rat IgG (STAR69...)

[FITC](#)

Goat Anti Rat IgG (STAR72...)

[HRP](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M315297:180503'

Printed on 05 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)