

Datasheet: MCA1553F

Description:	RAT ANTI HUMAN INTERLEUKIN-2:FITC
Specificity:	IL-2
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	MQ1-17H12
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry (1)	■			Neat

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) **Membrane permeabilization is required for this application. Bio-Rad recommend the use of Leucoperm™ (Product Code [BUF09](#)) for this purpose.**

Target Species	Human		
Species Cross Reactivity	Reacts with: Cynomolgus monkey N.B. Antibody reactivity and working conditions may vary between species.		
Product Form	Purified IgG - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant.		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml		
External Database Links	UniProt:		

Entrez Gene:

[3558](#) IL2 [Related reagents](#)

Specificity	Rat anti Human Interleukin-2 antibody, clone MQ1-17H12 recognizes human interleukin-2 (IL-2), also known as T-cell growth factor or Aldesleukin. IL-2 is a 153 amino acid secreted cytokine produced by T-cells following antigenic or mitogenic stimulation. It plays a crucial role in T-cell proliferation and immune regulation. Rat anti Human Interleukin-2 antibody, clone MQ1-17H12 has been successfully used for the detection and measurement of native and recombinant human IL-2 using flow cytometry and ELISA where it acts as the capture reagent in a very sensitive assay for IL-2 determination in human samples. IL-2 is used extensively for the treatment of patients with metastatic melanoma (Atkins et al. 1999) and renal cancer with reproducible beneficial effects (Rosenberg 2014). Monitoring of IL-2 levels in such patients is required for optimization of treatment regimens (Lotze et al. 1986). Rat anti Human Interleukin-2 antibody, clone MQ1-17H12 is also useful for evaluating levels of IL-2 in patients with autoimmune conditions including lupus erythematosus and rheumatoid arthritis where there may be an impaired ability to produce IL-2 (Talal et al. 1983).
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
References	1. Garbe, Y. <i>et al.</i> (2011) Semiallogenic fusions of MSI(+) tumor cells and activated B cells induce MSI-specific T cell responses. BMC Cancer. 11: 410. 2. Moncrieffe, H. <i>et al.</i> (2010) High expression of the ectonucleotidase CD39 on T cells from the inflamed site identifies two distinct populations, one regulatory and one memory T cell population. J Immunol. 185 (1): 134-43.
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL:FITC \(MCA6005F\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com 'M313633:180328'	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

Printed on 05 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)