

Datasheet: MCA193

Description:	MOUSE ANTI RAT IgE
Specificity:	IgE
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	MARE-1
Isotype:	IgG1
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			5ug/ml as coating antibody
Immunoprecipitation			■	
Western Blotting	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on IR2 IgE rat myeloma from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.1% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Rat IR162, IR1016, IR2 and IR410 IgE myeloma proteins.
External Database Links	UniProt: P01855 Related reagents

Entrez Gene:[299351](#) Ighe [Related reagents](#)

Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line.
Specificity	Mouse anti Rat IgE antibody, clone MARE-1 recognizes rat epsilon heavy chain of immunoglobulin IgE and does not cross-react with other classes of rat immunoglobulin. Mouse anti Rat IgE antibody, clone MARE-1 binds to rat IgE with an avidity of $4 \times 10^9 \text{M}^{-1}$
ELISA	This antibody may be used as a coating antibody in a sandwich ELISA in combination with clone MARK-1/MARL-15 (Product code MCA1296P) as detection reagent and purified rat IgE (Product code PRP07A).
References	1. Negrão-Corrêa, D. <i>et al.</i> (1996) Intestinal transport and catabolism of IgE: a major blood-independent pathway of IgE dissemination during a <i>Trichinella spiralis</i> infection of rats. J Immunol. 157 (9): 4037-44. 2. Bazin, H. <i>et al.</i> (1984) Rat monoclonal antibodies. I. Rapid purification from <i>in vitro</i> culture supernatants. J Immunol Methods. 66 (2): 261-9. 3. Bazin, H. <i>et al.</i> (1974) Transplantable immunoglobulin-secreting tumours in rats. IV. Sixty-three IgE-secreting immunocytoma tumours. Immunology. 26 (4): 713-23. 4. Bazin, H. <i>et al.</i> (1978) Transplantable IgD immunoglobulin-secreting tumors in rat. J Immunol. 121 (5): 2077-82. 5. Bazin, H. <i>et al.</i> (1984) Rat monoclonal antibodies. II. A rapid and efficient method of purification from ascitic fluid or serum. J Immunol Methods. 71 (1): 9-16. 6. Cho, J.K. & Cho, S.W. (2000) Shared epitope for monoclonal IR162 between <i>Anisakis simplex</i> larvae and <i>Clonorchis sinensis</i> and cross-reactivity between antigens. J Parasitol. 86 (5): 1145-9. 7. Silveira, M.R. <i>et al.</i> (2002) Infection with <i>Strongyloides venezuelensis</i> induces transient airway eosinophilic inflammation, an increase in immunoglobulin E, and hyperresponsiveness in rats. Infect Immun. 70: 6263-72. 8. Korinek, M. <i>et al.</i> (2016) Anti-allergic potential of <i>Typhonium blumei</i>: inhibition of degranulation via suppression of PI3K/PLCγ2 phosphorylation and calcium influx. Phytomedicine. 23 (14): 1706-15. 9. Bąbołowska, E. & Brzezińska-błaszczak, E. (2015) Human-derived cathelicidin LL-37 directly activates mast cells to proinflammatory mediator synthesis and migratory response. Cell Immunol. 293 (2): 67-73. 10. Agier, J. <i>et al.</i> (2018) Cathelicidin LL-37 Affects Surface and Intracellular Toll-Like Receptor Expression in Tissue Mast Cells. J Immunol Res. 2018: 7357162.
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10303 available at: 10303: https://www.bio-rad-antibodies.com/uploads/MSDS/10303.pdf

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG1 (HCA036...)	HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA1209\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M314178:180412'

Printed on 16 May 2018