

Datasheet: MCA2044F

Description:	MOUSE ANTI HUMAN HLA G:FITC
Specificity:	HLA G
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	MEM-G/9
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Preparation	Purified IgG prepared by affinity chromatography on Protein A		
Buffer Solution	Phosphate buffered saline		
Preservative Stabilisers	0.09% Sodium Azide		
	1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml		
Immunogen	Recombinant human HLA-G refolded with beta 2 microglobulin.		
External Database Links	UniProt:		
	P17693 Related reagents		
	Entrez Gene:		
	3135 HLA-G Related reagents		

Synonyms	HLA-6.0, HLAG
Fusion Partners	Spleen cells from immunised Balb/c mice were fused with myeloma cells.
Specificity	Mouse anti Human HLA G antibody, clone MEM-G/9 recognizes human HLA-G, a non-classical major histocompatibility complex (MHC) molecule. HLA-G expression is restricted to trophoblast cells and some medullary thymic epithelial cells. Several isoforms of the HLA-G molecule exist, which include the membrane bound isoforms HLA-G1 – G4 and soluble isoforms HLA-G5 – G7. Clone MEM-G/9 specifically recognizes surface expressed native HLA-G1, when associated with beta 2 microglobulin, but not does recognize the isoforms HLA-G2, G3 and G4. CMouse anti Human HLA G antibody, clone MEM-G/9 has also been reported to recognize the soluble isoform HLA-G5.
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	1. Fournel, S. <i>et al.</i> (2000) Comparative reactivity of different HLA-G monoclonal antibodies to soluble HLA-G molecules. Tissue Antigens. 55 (6): 510-8. 2. Menier, C. <i>et al.</i> (2003) Characterization of monoclonal antibodies recognizing HLA-G or HLA-E: new tools to analyze the expression of nonclassical HLA class I molecules. Hum Immunol. 64 (3): 315-26. 3. Kotze, D.J. <i>et al.</i> (2010) Embryo selection criteria based on morphology VERSUS the expression of a biochemical marker (sHLA-G) and a graduated embryo score: prediction of pregnancy outcome. J Assist Reprod Genet. 27 (6): 309-16. 4. Guetta, E. <i>et al</i> (2005) Trophoblasts isolated from the maternal circulation: <i>in vitro</i> expansion and potential application in non-invasive prenatal diagnosis. J Histochem Cytochem. 53: 337-9. 5. Hiby, S.E. <i>et al</i> (2010) Maternal activating KIRs protect against human reproductive failure mediated by fetal HLA-C2 J Clin Invest. 120: 4102-10. 6. Sher, G. <i>et al.</i> (2005) Influence of early ICSI-derived embryo sHLA-G expression on pregnancy and implantation rates: a prospective study. Hum Reprod. 20: 1359-63. 7. Sher, G. <i>et al.</i> (2005) Soluble human leukocyte antigen G expression in phase I culture media at 46 hours after fertilization predicts pregnancy and implantation from day 3 embryo transfer. Fertil Steril. 83: 1410-3. 8. Apps, R. <i>et al.</i> (2011) Ex vivo functional responses to HLA-G differ between blood and decidual NK cells. Mol Hum Reprod. 17: 577-86. 9. Manaster, I. <i>et al.</i> (2012) MiRNA-mediated control of HLA-G expression and function. PLoS One. 7: e33395. 10. Nüchel, H. <i>et al.</i> (2005) HLA-G expression is associated with an unfavorable outcome and immunodeficiency in chronic lymphocytic leukemia. Blood. 105: 1694-8. 11. Yao, Y.Q. <i>et al.</i> (2005) Differential expression of alternatively spliced transcripts of HLA-G in human preimplantation embryos and inner cell masses. J Immunol. 175 (12): 8379-85. 12. de Carvalho, J.F. <i>et al.</i> (2012) Heparin increases HLA-G levels in primary antiphospholipid syndrome. Clin Dev Immunol. 2012: 232390. 13. Guenther, S. <i>et al.</i> (2012) Decidual macrophages are significantly increased in spontaneous miscarriages and over-express FasL: a potential role for macrophages in trophoblast apoptosis. Int J Mol Sci. 13 (7): 9069-80. 14. Apps, R. <i>et al.</i> (2011) Genome-wide expression profile of first trimester villous and extravillous human trophoblast cells. Placenta. 32 (1): 33-43. 15. Lim DS <i>et al.</i> (2014) The combination of type I IFN, TNF-α, and cell surface receptor engagement with dendritic cells enables NK cells to overcome immune evasion by dengue virus. J Immunol. 193 (10): 5065-75. 16. Reches, A. <i>et al.</i> (2016) HNRNPR Regulates the Expression of Classical and Nonclassical MHC Class I Proteins. J Immunol. 196 (12): 4967-76.

17. Bröker, P *et al.* (2012) A nanostructured SAW chip-based biosensor detecting cancer cells
[Sensors and Actuators B: Chemical. 165 \(1\): 1-6.](#)

Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:FITC \(MCA928F\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M301902:170109'

Printed on 01 May 2018