

Datasheet: MCA2164PE

Description:	MOUSE ANTI CHICKEN CD4:RPE
Specificity:	CD4
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	2-35
Isotype:	IgG2b
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat - 1/5

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Chicken								
Species Cross Reactivity	Reacts with: Turkey N.B. Antibody reactivity and working conditions may vary between species.								
Product Form	Purified IgG - lyophilized								
Reconstitution	Reconstitute with 1.0 ml distilled water								
Max Ex/Em	<table border="1"> <thead> <tr> <th>Fluorophore</th> <th>Excitation Max (nm)</th> <th>Emission Max (nm)</th> </tr> </thead> <tbody> <tr> <td>RPE 488nm laser</td> <td>496</td> <td>578</td> </tr> </tbody> </table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	RPE 488nm laser	496	578		
Fluorophore	Excitation Max (nm)	Emission Max (nm)							
RPE 488nm laser	496	578							
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant								
Buffer Solution	Phosphate buffered saline								
Preservative Stabilisers	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose								
Immunogen	Chicken embryonic thymocytes								
Fusion Partners	Lymph node cells from immunised Balb/c mice were fused with cells of the SP2/0 myeloma cell line								

Specificity	Mouse anti Chicken CD4, clone 2-35 recognizes the chicken homologue of human CD4, a ~64 kDa cell surface protein expressed by thymocytes and a subset of T cells (Luhtala et al. 1993). Mouse anti Chicken CD4, clone 2-35 has been demonstrated to recognise turkey CD4 (Li et al. 1998).
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	<ol style="list-style-type: none"> 1. Vainio, O. <i>et al.</i> (1989) Characterization of chicken CD4-expressing cells. Prog Clin Biol Res. 307: 45-56. 2. Li, Z. <i>et al.</i> (1999) Cross-reactive anti-chicken CD4 and CD8 monoclonal antibodies suggest polymorphism of the turkey CD8alpha molecule. Poult Sci. 78 (11): 1526-31. 3. Koskinen, R. <i>et al.</i> (1999) Cloning and modeling of the first nonmammalian CD4. J Immunol. 162 (7): 4115-21. 4. Pavlova, S.P. <i>et al.</i> (2010) <i>In vitro</i> and <i>in vivo</i> characterization of glycoprotein C-deleted infectious laryngotracheitis virus. J Gen Virol. 91 (Pt 4): 847-57. 5. Luhtala, M. <i>et al.</i> (1993) Analysis of chicken CD4 by monoclonal antibodies indicates evolutionary conservation between avian and mammalian species. Hybridoma. 12: 633-46. 6. Rosa, A.C. <i>et al.</i> (2014) Isolation and molecular characterization of Brazilian turkey reovirus from immunosuppressed young poult. Arch Virol. 159 (6): 1453-7. 7. Blohm, U. <i>et al.</i> (2016) Immunological Competence of Different Domestic Chicken Breeds Against Avian Influenza Infection. Avian Dis. 60 (1 Suppl): 262-8. 8. Röhe, I. <i>et al.</i> (2017) Effect of feeding soybean meal and differently processed peas on the gut mucosal immune system of broilers Poultry Science. Feb 23 [Epub ahead of print] 9. Abd El-Hack, M. & Alagawany, M. (2015) Performance, egg quality, blood profile, immune function, and antioxidant enzyme activities in laying hens fed diets with thyme powder Journal of Animal and Feed Sciences. 24 (2): 127-33. 10. Kannan, T.A. <i>et al.</i> (2017) Age Related Changes in T Cell Subsets in Thymus and Spleen of Layer Chicken (<i>Gallus domesticus</i>) Intl J Curr Microbiol App Sci. 6 (1): 15-9. 11. Sachan, S. <i>et al.</i> (2015) Adjuvant potential of resiquimod with inactivated Newcastle disease vaccine and its mechanism of action in chicken. Vaccine. 33 (36): 4526-32.
Storage	<p>Prior to reconstitution store at +4°C. Following reconstitution store at +4°C.</p> <p>DO NOT FREEZE.</p> <p>This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10075 available at: 10075: https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL:RPE \(MCA691PE\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Printed on 05 May 2018
