

Datasheet: MCA2155

Description:	MOUSE ANTI HUMAN CD206
Specificity:	CD206
Other names:	MANNOSE RECEPTOR C TYPE 1
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	15-2
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen	▪			
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting (1)	▪			
Immunofluorescence	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1)**Clone 15-2 recognises a protein of approximately 175kDa under non-reducing conditions.**

Target Species	Human
Product Form	Purified IgG - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Purified human mannose receptor.
External Database Links	UniProt: P22897 Related reagents

Entrez Gene:[4360](#) MRC1 [Related reagents](#)**Synonyms**

CLEC13D

Fusion Partners

Spleen cells from immunised Balb/c mice where fused with cells of the SP2/0 Ag.14 mouse myeloma cell line.

Specificity

Mouse anti Human CD206 monoclonal antibody, clone 15-2 recognizes human macrophage mannose receptor 1, also known as CD206 or C-type lectin domain family 13 member D-like. CD206 is a ~175 kDa single pass type I transmembrane glycoprotein belonging to the group of pattern recognition receptors ([Paveley et al. 2011](#)). CD206 has multiple carbohydrate recognition motifs and acts as a receptor for bacteria, fungi and other pathogens ([Ezekowitz et al. 1990](#)). CD206 is predominantly expressed in tissue macrophages and dendritic cells ([Engering et al. 1997](#)) and can also be found in a subpopulation of endothelial cells ([Pack et al. 2007](#)) and sperm cells ([Cardona-Maya et al. 2006](#)). CD206 can also be detected in a soluble form in human plasma and is elevated in patients with acute sepsis ([Rødgaard-Hansen et al. 2013](#)).

Mouse anti CD206, clone 15-2 has been used extensively to monitor mannose receptor modulation in macrophages treated with a wide range of cytokines and growth factors ([Chang et al. 2004](#)) and to indicate CD206 as a marker for alternative activation of macrophages ([Joerink et al. 2011](#)).

Flow Cytometry

Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

1. Barrett-Bergshoeff, M. *et al.* (1997) Monoclonal antibodies against the human mannose receptor that inhibit the binding of tissue-type plasminogen activator. [Thromb Haemost. 77: 718-24.](#)
2. Koning, N. *et al.* (2009) Distribution of the immune inhibitory molecules CD200 and CD200R in the normal central nervous system and multiple sclerosis lesions suggests neuron-glia and glia-glia interactions. [J Neuropathol Exp Neurol. 68: 159-67.](#)
3. Emara, M. *et al.* (2011) Recognition of the major cat allergen Fel d 1 through the cysteine-rich domain of the mannose receptor determines its allergenicity. [J Biol Chem. 286:13033-40.](#)
4. Chang, S.K. *et al.* (2008) B lymphocyte stimulator regulates adaptive immune responses by directly promoting dendritic cell maturation. [J Immunol. 180: 7394-403.](#)
5. MacKinnon, A.C. *et al.* (2008) Regulation of alternative macrophage activation by galectin-3. [J Immunol. 180: 2650-8.](#)
6. Lai, W.K. *et al.* (2006) Expression of DC-SIGN and DC-SIGNR on human sinusoidal endothelium: a role for capturing hepatitis C virus particles. [Am J Pathol. 169: 200-8.](#)
7. Ueno, N. *et al.* (2009) Differences in human macrophage receptor usage, lysosomal fusion kinetics and survival between logarithmic and metacyclic *Leishmania infantum chagasi* promastigotes. [Cell Microbiol. 11: 1827-41.](#)
8. Kato, M. *et al.* (2000) Expression of multilectin receptors and comparative FITC-dextran uptake by human dendritic cells. [Int. Immunol. 12:1511-9.](#)
9. Chang, Y.C. *et al.* (2004) Modulation of macrophage differentiation and activation by decoy receptor 3. [J Leukoc Biol. 75: 486-94.](#)
10. Yamamoto, H. *et al.* (2011) Sphingosylphosphorylcholine and lysosulfatide have inverse regulatory functions in monocytic cell differentiation into macrophages. [Arch Biochem Biophys. 506: 83-91.](#)
11. He, L.Z. *et al.* (2007) Antigenic targeting of the human mannose receptor induces tumor immunity. [J Immunol. 178: 6259-67.](#)
12. Sturge, J. *et al.* (2007) Mannose receptor regulation of macrophage cell migration. [J Leukoc Biol. 82: 585-93.](#)
13. Torrelles, J.B. *et al.* (2006) Fine discrimination in the recognition of individual species of

- phosphatidyl-myo-inositol mannosides from *Mycobacterium tuberculosis* by C-type lectin pattern recognition receptors. [J Immunol. 177 \(3\): 1805-16.](#)
14. Larsson, K. *et al.* (2015) COX/mPGES-1/PGE2 pathway depicts an inflammatory-dependent high-risk neuroblastoma subset. [Proc Natl Acad Sci U S A. 112 \(26\): 8070-5.](#)
15. Stankevich, K.S. *et al.* (2015) Surface modification of biomaterials based on high-molecular polylactic acid and their effect on inflammatory reactions of primary human monocyte-derived macrophages: perspective for personalized therapy. [Mater Sci Eng C Mater Biol Appl. 51: 117-26.](#)
16. Yamane, K. & Leung, K.P. (2016) Rabbit M1 and M2 macrophages can be induced by human recombinant GM-CSF and M-CSF. [FEBS Open Bio. 6 \(9\): 945-53.](#)
17. Karna, S.L. *et al.* (2016) RNA-Seq Transcriptomic Responses of Full-Thickness Dermal Excision Wounds to *Pseudomonas aeruginosa* Acute and Biofilm Infection. [PLoS One. 11 \(10\): e0165312.](#)
18. Özçelik H *et al.* (2015) Harnessing the multifunctionality in nature: a bioactive agent release system with self-antimicrobial and immunomodulatory properties. [Adv Healthc Mater. 4 \(13\): 2026-36.](#)
19. Argueta-Donohué J *et al.* (2016) Differential phagocytosis of *Leishmania mexicana* promastigotes and amastigotes by monocyte-derived dendritic cells. [Microbiol Immunol. 60 \(6\): 369-81.](#)

Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG1 (HCA036...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL (MCA928)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M318515:180718'

Printed on 01 Aug 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)