

Datasheet: MCA2144EL

Description:	MOUSE ANTI HUMAN CD253:Low Endotoxin
Specificity:	CD253
Other names:	TRAIL
Format:	Low Endotoxin
Product Type:	Monoclonal Antibody
Clone:	2E5
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			
Functional Assays	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Endotoxin Level	<0.01EU/ug
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Recombinant human CD253 (TRAIL) aa 95-281.

**External Database
Links**

UniProt:

[P50591](#)

[Related reagents](#)

Entrez Gene:

[8743](#)

TNFSF10

[Related reagents](#)

Synonyms

APO2L, TRAIL

Specificity

Mouse anti Human CD253 antibody, clone 2E5 recognizes the 281 amino acid protein termed TNF related apoptosis-inducing ligand (TRAIL), also known as Apo2L.

TRAIL is a cytotoxic protein closely related to Fas (CD95), which induces extensive apoptosis in a range of cell types. TRAIL induces apoptosis via the receptors DR4 and DR5. Most normal tissues are resistant to TRAIL induced apoptosis due to their expression of decoy receptors, DcR1 and DcR2, which antagonise TRAIL activity.

TRAIL is expressed by cell lines, activated lymphocytes and on syncytiotrophoblast cells in first trimester pregnancy.

This antibody neutralizes TRAIL activity.

References

1. Plasilova, M. *et al.* (2002) TRAIL (Apo2L) suppresses growth of primary human leukemia and myelodysplasia progenitors. [Leukemia. 16: 67-73.](#)

Storage

Store at -20°C only.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

**Health And Safety
Information**

Material Safety Datasheet Documentation #10162 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP

Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Human Anti Mouse IgG1 (HCA036...)	HRP
Goat Anti Mouse IgG (STAR70...)	FITC

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M280828:160205'

Printed on 07 Jun 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)