

Datasheet: MCA2307F

Description:	MOUSE ANTI PIG CD5:FITC
Specificity:	CD5
Format:	FITC
Product Type:	Monoclonal Antibody
Clone:	1H6/8
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10
Immunoprecipitation	▪			
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Pig		
Product Form	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	FITC	490	525
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative Stabilisers	0.09% Sodium Azide		
	1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.1 mg/ml		
Immunogen	ConA/PMA activated porcine peripheral blood cells.		
External Database Links	UniProt:		
	Q9GMA0	Related reagents	

Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse X63-Ag.8.653 myeloma cell line.
Specificity	Mouse anti Pig CD5, clone 1H6/8 recognizes pig CD5. Clone 1H6/8 was clustered as CD5a at the Second International Swine CD antigen Workshop (Saalmuller et al. 1998), based on its cellular distribution, ability to completely inhibit the binding of monoclonal antibody clone b53b7, the internal control for the second workshop and the immunoprecipitation of the same 55 kDa + 60 kDa doublet as the control antibody (Pescowitz et al. 1998). Porcine CD5 is expressed by all circulating peripheral blood T lymphocytes and some B cells. CD5 is not expressed by porcine $\gamma\delta$ T cells (Saalmuller et al. 1994) and can hence be used for the discrimination between NK cells (CD4- CD8+ CD5-) and MHC-restricted cytotoxic T cells (CD4- CD8+ CD5+). No reactivity of Mouse anti pig CD5 was noted with canine, bovine, equine or human tissue samples (Doménech et al. 2003).
Flow Cytometry	Use 10ul of the suggested working dilution to 1×10^6 cells in 100ul.
References	1. Pescovitz, M.D. <i>et al.</i> (1998) Analyses of monoclonal antibodies reacting with porcine CD5: results from the Second International Swine CD Workshop. Vet Immunol Immunopathol. 60 (3-4): 269-73. 2. Saalmüller, A. <i>et al.</i> (1998) Overview of the Second International Workshop to define swine cluster of differentiation (CD) antigens. Vet Immunol Immunopathol. 60 (3-4): 207-28. 3. Doménech, N. <i>et al.</i> (2003) A new epitope on swine CD5 molecule detected by monoclonal antibody 5F12/9. Hybrid Hybridomics. 22 (3): 179-82. 4. Sánchez-torres, C. <i>et al.</i> (2003) Expression of porcine CD163 on monocytes/macrophages correlates with permissiveness to African swine fever infection. Arch Virol. 148 (12): 2307-23.
Further Reading	1. Piriou-Guzylack, L. (2008) Membrane markers of the immune cells in swine: an update. Vet Res. 39: 54.
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost free freezers is not recommended. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL:FITC \(MCA929F\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com "M301689:170109"	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

Printed on 05 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)