

Datasheet: MCA2421PE

Description:	MOUSE ANTI HUMAN CD244:RPE
Specificity:	CD244
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	2B4.69
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
Reconstitution	Reconstitute with 1.0 ml distilled water		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1% Bovine Serum Albumin		
	5% Sucrose		
Immunogen	CD244 Fusion protein.		

External Database Links

UniProt:

[Q9BZW8](#)

[Related reagents](#)

Entrez Gene:

Synonyms	2B4
-----------------	-----

Fusion Partners	Spleen cells from immunised Balb/c mice were fused with cells of the NS-1 myeloma cell line.
------------------------	--

Specificity	<p>Mouse anti Human CD244 antibody, clone 2B4.69 recognizes the CD244 cell surface antigen, also known as natural killer receptor 2B4 (NKG2B4). CD244 is a 370 amino acid single pass type I transmembrane glycoprotein with two Ig-like V type domains. CD244 is a member of the SLAM family of receptors which is expressed on natural killer (NK) cells, gamma delta T cells, CD14+ monocytes, and a subset of CD8+ T cells.</p> <p>CD244 was originally identified as an activating receptor on NK cells, although recent studies suggest that the receptor may have a wide range of activities including NK cell activation, costimulation, or even inhibition, depending on the cell type and stage of differentiation. The ligand of CD244 is CD48, a GPI-anchored membrane protein also belonging to the SLAM family.</p>
--------------------	--

Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
-----------------------	---

References	<ol style="list-style-type: none">1. Romero, X. <i>et al.</i> (2005) CD229 (Ly9) lymphocyte cell surface receptor interacts homophilically through its N-terminal domain and relocalizes to the immunological synapse. J Immunol. 174 (11): 7033-42.2. Bengsch, B. <i>et al.</i> (2014) Restoration of HBV-specific CD8+ T cell function by PD-1 blockade in inactive carrier patients is linked to T cell differentiation. J Hepatol. 61 (6): 1212-9.3. Raziorrouh, B. <i>et al.</i> (2010) The Immunoregulatory Role of CD244 in Chronic Hepatitis B Infection and its Inhibitory Potential on Virus-Specific CD81 T-cell Function. Hepatology. 52: 1934-47.4. Boussaad, I. <i>et al.</i> (2011) Wild-type measles virus interferes with short-term engraftment of human CD34+ hematopoietic progenitor cells. J Virol. 85: 7710-8.
-------------------	--

Storage	<p>Prior to reconstitution store at +4°C. After reconstitution store at +4°C. DO NOT FREEZE.</p> <p>This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	--

Shelf Life	18 months from date of reconstitution.
-------------------	--

Health And Safety Information	Material Safety Datasheet documentation #10075 available at: 10075: https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf
--------------------------------------	---

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA928PE\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com
'M300761:170106'

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

Printed on 05 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)