

Datasheet: MCA2464PBT

Description:	RAT ANTI MOUSE CD279:Pacific Blue®
Specificity:	CD279
Other names:	PD-1
Format:	Pacific Blue®
Product Type:	Monoclonal Antibody
Clone:	RMP1-30
Isotype:	IgG2b
Quantity:	25 TESTS/0.25ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/5

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse						
Product Form	Purified IgGconjugated to Pacific Blue® - liquid						
Max Ex/Em	<table border="1"> <thead> <tr> <th>Fluorophore</th> <th>Excitation (nm)</th> <th>Emission (nm)</th> </tr> </thead> <tbody> <tr> <td>Pacific Blue®</td> <td>410</td> <td>455</td> </tr> </tbody> </table>	Fluorophore	Excitation (nm)	Emission (nm)	Pacific Blue®	410	455
Fluorophore	Excitation (nm)	Emission (nm)					
Pacific Blue®	410	455					
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant						
Buffer Solution	Phosphate buffered saline						
Preservative	0.09% Sodium Azide						
Stabilisers	1% Bovine Serum Albumin						
Approx. Protein Concentrations	IgG concentration 0.05 mg/ml						
Immunogen	PD-1 transferred BHK cells.						
External Database Links	UniProt: Q02242 Related reagents						

Entrez Gene:

[18566](#) Pcd1 [Related reagents](#)

Synonyms	Pd1
-----------------	-----

Fusion Partners	Spleen cells from immunised Sprague Dawley rats were fused with cells of the P3U1 myeloma cell line.
------------------------	--

Specificity	<p>Rat anti Mouse CD279 antibody, clone RMP1-30 recognizes mouse CD279, a ~55 kDa cell surface protein, a member of the CD28/CTLA-4 family, otherwise known as Programmed Death-1 (PD-1). CD279 is expressed predominantly on activated T- and B- lymphocytes and on a subset of thymocytes.</p> <p>Studies suggest that CD279, an immunoinhibitory receptor, plays a critical role in peripheral tolerance induction and prevention of autoimmune disease. Two members of the B7 family, B7-H1 (PD-L1) and B7-DC (PD-L2), have been identified as the ligands for CD279.</p> <p>Rat anti Mouse CD279 antibody, clone RMP1-30 does not block the binding of either B7-H1-Ig or B7-DC-Ig fusion proteins to PD-1 transfected BHK cells.</p>
--------------------	---

Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
-----------------------	---

References	<ol style="list-style-type: none">1. Matsumoto, K. <i>et al.</i> (2004) B7-DC regulates asthmatic response by an IFN-gamma-dependent mechanism. J Immunol. 172 (4): 2530-41.2. Bartkowiak, T. (2013) Novel Imaging-Based Techniques Reveal a Role for PD-1/PD-L1 in Tumor Immune Surveillance in the Lung. UT GSBS Dissertations and Theses (Open Access). Paper 354.3. Hu, Z. <i>et al.</i> (2013) Regulatory CD8+ T cells associated with erosion of immune surveillance in persistent virus infection suppress <i>in vitro</i> and have a reversible proliferative defect. J Immunol. 191 (1): 312-22.
-------------------	---

Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. This product is photosensitive and should be protected from light.</p> <p>Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	--

Shelf Life	18 months from date of despatch.
-------------------	----------------------------------

Acknowledgements	This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchased product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com
-------------------------	---

Health And Safety	Material Safety Datasheet Documentation #10041 available at:
--------------------------	--

Information <https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf>

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE SEROBLOCK FcR \(BUF041A\)](#)

[MOUSE SEROBLOCK FcR \(BUF041B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

 Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

 Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

 Email: antibody_sales_de@bio-rad.com

'M287843:160607'

Printed on 07 Jun 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)