

Datasheet: MCA2472

Description:	MOUSE ANTI PHOSPHOTYROSINE
Specificity:	PHOSPHOTYROSINE
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	PY20
Isotype:	IgG2b
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			100ng/ml
Immunoprecipitation			■	
Western Blotting (1)	■			2ug/ml
Immunofluorescence	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) **For the detection of phosphoproteins, tyrosine phosphatase inhibitors such as 1mM Sodium Orthovanadate should be added to the sample buffer. Milk or other casein-based blocking solutions are not recommended as casein is a phosphoprotein and its use can result in high background.**

Target Species	Chemical
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography.
Buffer Solution	Borate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Specificity	Mouse anti Phosphotyrosine antibody, clone PY20 recognizes phosphotyrosine, enabling the

detection, characterisation and isolation of proteins containing phosphorylated tyrosine residues.

The phosphorylation of tyrosine acts as a important signal in the control of cell mitogenesis, differentiation, proliferation, and migration and occurs following the activation of intracellular tyrosine kinases, including the T-cell receptor (TCR), epidermal growth factor (EGF) and many families of receptor and non-receptor protein tyrosine kinases (PTKs), which catalyse the transfer of ATP to a tyrosine residue on specific cell protein targets.

The binding of PY20 to phosphorylated tyrosines can be inhibited by free phosphotyrosine and phenylphosphate, but not by free phosphate, phosphoserine or phosphothreonine.

The affinity of PY20 for phosphotyrosine is 10^{-6} to 10^{-7} M.

References

1. Ruff-Jamison, S. *et al.* (1991) Heavy and light chain variable region sequences and antibody properties of anti-phosphotyrosine antibodies reveal both common and distinct features. [J Biol Chem. 266 \(10\): 6607-13.](#)
2. Takagi, S. *et al.* (1991) Intracellular localization of tyrosine kinase substrates beneath crosslinked surface immunoglobulins in B cells. [J Exp Med. 174 \(2\): 381-8.](#)
3. Glenney, J.R. Jr. *et al.* (1988) Monoclonal antibodies to phosphotyrosine. [J Immunol Methods. 109 \(2\): 277-85.](#)
4. Vendel, A.C. *et al.* (2009) B and T lymphocyte attenuator regulates B cell receptor signaling by targeting Syk and BLNK. [J Immunol. 182: 1509-17.](#)

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

Health And Safety Information

Material Safety Datasheet documentation available at:
Material Safety Datasheet Documentation #10077 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10077.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR76...)	RPE
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Human Anti Mouse IgG2b (HCA038...)	FITC , HRP
Goat Anti Mouse IgG (STAR70...)	FITC
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP

Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M251140:140326'

Printed on 04 Oct 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)