

Datasheet: MCA2514A647T

Description:	MOUSE ANTI HUMAN CD47:Alexa Fluor® 647
Specificity:	CD47
Other names:	SIRP LIGAND
Format:	ALEXA FLUOR® 647
Product Type:	Monoclonal Antibody
Clone:	1/1A4
Isotype:	IgG2b
Quantity:	25 TESTS/0.25ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat - 1/5

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to Alexa Fluor® 647 - liquid		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	Alexa Fluor®647	650	665
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide (NaN ₃)		
Stabilisers	1% Bovine Serum Albumin		
Approx. Protein Concentrations	IgG concentration 0.05mg/ml		
Immunogen	Human erythrocytes.		
External Database Links	UniProt:		
	Q08722	Related reagents	

Entrez Gene:

[961](#) CD47 [Related reagents](#)

Synonyms	MER6
-----------------	------

Fusion Partners	Spleen cells from immunised Balb/c mice were fused with cells of the mouse NS-1 myeloma cell line.
------------------------	--

Specificity	<p>Mouse anti Human CD47 antibody, clone 1/1A4 recognizes human CD47, also known as integrin-associated protein, a 50kDa glycoprotein expressed on all haematopoietic cell lines, epithelial cells, endothelial cells, fibroblasts and many tumour cell lines.</p> <p>CD47 is involved in several cellular processes including regulation of T cell and dendritic cell activation, regulation of neutrophil (PMN) chemotaxis, caspase-independent apoptosis, and may act as a signal transducer in the regulation of cation fluxes across cell membranes. CD47 has been identified as a ligand for the inhibitory macrophage receptor known as signal regulatory protein alpha (SIRPalpha) and as a receptor for thrombospondin-1.</p> <p>Cross-linking of CD47 with clone 1/1A4 co-stimulates the proliferation of T cells (Waclavicek et al. 1993).</p>
--------------------	--

Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
-----------------------	---

References	1. Waclavicek, M. <i>et al.</i> (1997) T cell stimulation via CD47: agonistic and antagonistic effects of CD47 monoclonal antibody 1/1A4. J Immunol. 159 (11): 5345-54.
-------------------	---

Further Reading	1. Frazier, W.A. <i>et al.</i> (1999) The thrombospondin receptor integrin-associated protein (CD47) functionally couples to heterotrimeric Gi. J Biol Chem. 274 (13): 8554-60.
------------------------	---

Storage	<p>Store at +4°C or at -20°C if preferred.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product is photosensitive and should be protected from light.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
----------------	--

Shelf Life	18 months from date of despatch.
-------------------	----------------------------------

Acknowledgements	<p>This product is provided under an intellectual property licence from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchase product solely in research, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad CA 92008 USA or outlicensing@thermofisher.com</p>
-------------------------	---

Health And Safety Information	Material Safety Datasheet documentation #10041 available at: 10041: https://www.bio-rad-antibodies.com/uploads/MSDS/10041.pdf
--------------------------------------	---

Related Products

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL:Alexa Fluor® 647 \(MCA691A647\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M312206:180209'

Printed on 20 Jun 2018