

Datasheet: MCA2586PE

Description:	RAT ANTI MOUSE CDw199:RPE
Specificity:	CDw199
Other names:	CCR9
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	9B1
Isotype:	IgG2a
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse		
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
Reconstitution	Reconstitute with 1.0ml distilled water		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578
Preparation	Purified IgG prepared by affinity chromatography on Protein G from		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide (NaN ₃)		
Stabilisers	1%	Bovine Serum Albumin	
	5%	Sucrose	
Immunogen	25 amino acids from the amino terminal of murine CCR9, coupled to		
External Database Links	UniProt: Q9WUT7 Related reagents		

Entrez Gene:[12769](#) Ccr9 [Related reagents](#)

Synonyms	Cmkbr10
Fusion Partners	Spleen cells from immunised Lewis rat were fused with cells of the X63-Ag8.653 myeloma cell line.
Specificity	Rat anti Mouse CDw199 antibody, clone 9B1 detects mouse chemokine receptor 9 (CCR9), also known as CDw199. CCR9 is a seven transmembrane protein which is a member of the beta chemokines receptor family. CCR9 is expressed in thymus and small intestine, and mediates the recruitment and development of thymocytes in response to CCL25; the only known ligand of CCR9. CCR9 is an alternative co-receptor with CD4 for HIV infection, and has been implicated in adult T cell leukaemia. Two alternatively spliced transcript variants have been described. Clone 9B1 blocks CCL25-induced Ca²⁺ release, attenuating inflammation in a murine model of chronic ileitis when given in early disease stages.
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10⁶ cells in 100ul. The Fc region of monoclonal antibodies may bind non-specifically to cells expressing low affinity Fc receptors. This may be reduced by using SeroBlock FcR (BUF041A/B).
References	1. Rivera-nieves, J. <i>et al.</i> (2006) Antibody blockade of CCL25/CCR9 ameliorates early but not late chronic murine ileitis. Gastroenterology. 131 (5): 1518-29. 2. Nagakubo, D. <i>et al.</i> (2007) Expression of CCR9 in HTLV-1+ T cells and ATL cells expressing Tax. Int J Cancer. 120 (7): 1591-7.
Further Reading	1. Uehara, S. <i>et al.</i> (2006) Premature expression of chemokine receptor CCR9 impairs T cell development. J Immunol. 176 (1): 75-84.
Storage	Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of reconstitution
Health And Safety Information	Material Safety Datasheet documentation #10075 available at: 10075: https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

[RAT IgG2a NEGATIVE CONTROL:RPE \(MCA1212PE\)](#)

Recommended Useful Reagents

[MOUSE SEROBLOCK FcR \(BUF041A\)](#)
[MOUSE SEROBLOCK FcR \(BUF041B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M304026:170314'

Printed on 05 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)