

Datasheet: MCA2673GA

Description:	MOUSE ANTI AMPK ALPHA 1/2
Specificity:	AMPK ALPHA 1/2
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	34.2
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA	■			
Immunoprecipitation	■			
Western Blotting	■			1/100 - 1/1000
Immunofluorescence	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Drosophila
Species Cross Reactivity	Reacts with: Vertebrates, Human, Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml

Immunogen 29 amino acid peptide within the N-terminal catalytic domain of Drosophila AMPK alpha.

**External Database
Links**

UniProt:

[P54646](#) [Related reagents](#)
[Q13131](#) [Related reagents](#)

Entrez Gene:

[5563](#) PRKAA2 [Related reagents](#)
[5562](#) PRKAA1 [Related reagents](#)

Synonyms

AMPK, AMPK1, AMPK2

Specificity

Mouse anti AMPK Alpha 1/2 antibody, clone 34.2 recognises both the alpha 1 and the alpha 2 isoforms of AMP-activated protein kinase (AMPK). AMPK is a highly conserved heterotrimeric protein kinase found in all eukaryotes. The heterotrimeric complex exists as multiple isoforms and consists of a catalytic alpha subunit, which is stabilised by non-catalytic regulatory beta- and gamma- subunits. AMPK plays a vital role in maintaining cellular energy homeostasis and is activated by rising levels of AMP under conditions of hypoxia, ischemia, exercise, low glucose levels.

AMPK may play an important role in the control of type 2 diabetes and metabolic syndromes. Recent studies suggest that AMPK functions as a suppressor of cell proliferation of both malignant and normal cells.

Western Blotting

MCA2673GA detects a band of approximately 63kDa in rat brain cell lysates.

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE

Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR8...) [DyLight®800](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Human Anti Mouse IgG2a (HCA037...) [FITC](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M315829:180503'

Printed on 05 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)