

Datasheet: MCA2604

Description:	MOUSE ANTI HUMAN C3a
Specificity:	C3a
Other names:	COMPLEMENT COMPONENT 3a
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	10-21
Isotype:	IgG1
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			1/1500 - 1/6000
Western Blotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Borate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	1.0 mg/ml
Immunogen	Native purified human C3a.

External Database Links

UniProt:

[P01024](#) [Related reagents](#)

Entrez Gene:

[718](#) C3 [Related reagents](#)

Synonyms	CPAMD1
Specificity	Mouse anti Human C3a antibody, clone 10-21 recognizes the ~9 kDa complement component 3a (C3a) which is present in the blood serum. C3a is generated over the course of complement activation where complement 3 (C3) convertases cleave C3 to C3a and C3b, which is further degraded into iC3b and C3dg/C3d. C3 is crucial in the induction of tolerance generated when an antigen is introduced into immunoprivileged sites and this is exploited by pathogens and cancer cells to evade the immune system by inhibiting complement activation. C3a mediates local inflammatory process. It induces the contraction of smooth muscle, increases vascular permeability and causes histamine release from mast cells and basophilic leukocytes. C3a is expressed at significantly higher levels in colorectal cancer. Mouse anti Human C3a antibody, clone 10-21 shows minimal reactivity with other complement proteins in radioimmunoassays and haemagglutination. Mouse anti Human C3a antibody, clone 10-21 has been reported to work in western blotting applications.
References	1. Hugli, T.E. (1975) Human anaphylatoxin (C3a) from the third component of complement. Primary structure. J Biol Chem. 250 (21): 8293-301. 2. Habermann, J.K. <i>et al.</i> (2006) Increased serum levels of complement C3a anaphylatoxin indicate the presence of colorectal tumors. Gastroenterology. 131 (4): 1020-9; quiz 1284.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10334 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10334.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP

Rabbit Anti Mouse IgG (STAR8...) [DyLight®800](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Human Anti Mouse IgG1 (HCA036...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M311789:180104'

Printed on 04 Jan 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)