

Datasheet: MCA2666

Description:	MOUSE ANTI HUMAN C6
Specificity:	C6
Other names:	COMPLEMENT COMPONENT 6
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	056B-214.2.4.2
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA	■			
Immunoprecipitation			■	
Western Blotting			■	
Functional Assays			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Borate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Purified human C6.
External Database	UniProt:

Links

[P13671](#) [Related reagents](#)

Entrez Gene:

[729](#) C6 [Related reagents](#)

Specificity	Mouse anti Human C6 antibody, clone 056B-214.2.4.2 recognises complement component 6 (C6), a 104 kDa member of the complement C6/C7/C8/C9 family present in blood serum. The factor I modules of C6 bind the C-terminus of the C5b alpha-chain, anchoring the molecule to the membrane; they also binds complement component 7. This enables the formation of the C5b-7 precursor complex. C5b-7 then functions as a receptor for components C8 and C9 which form the complement membrane attack complex (MAC). Deficiency of C6 results in susceptibility to bacterial infection. Mouse anti Human C6 antibody, clone 056B-214.2.4.2 recognises C6 as a component of the MAC.
Histology Positive Control Tissue	Kidney from post streptococcal glomerulonephritis patients
References	1. DiScipio, R.G. <i>et al.</i> (1999) Function of the factor I modules (FIMS) of human complement component C6. J Biol Chem. 274 (45): 31811-8.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10334 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10334.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP

Human Anti Mouse IgG1 (HCA036...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Recommended Useful Reagents

[HISTAR DETECTION SYSTEM \(STAR3000A\)](#)

[HISTAR DETECTION SYSTEM \(STAR3000B\)](#)

[HISTAR DETECTION SYSTEM \(STAR3000C\)](#)

**North & South
America**

Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M312066:180117'

Printed on 19 Feb 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)