

Datasheet: MCA2645

Description:	MOUSE ANTI HUMAN C8
Specificity:	C8
Other names:	COMPLEMENT COMPONENT 8
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	056B-373
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA	■			
Immunoprecipitation			■	
Western Blotting	■			
Functional Assays			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A
Buffer Solution	Borate buffered saline
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Purified human C8
External Database	UniProt:

Links	P07360	Related reagents
	P07358	Related reagents
	P07357	Related reagents

Entrez Gene:

733	C8G	Related reagents
732	C8B	Related reagents
731	C8A	Related reagents

Specificity	Mouse anti Human C8 antibody, clone 056B-373 recognizes complement component 8 (C8), a 151 kDa member of the complement C6/C7/C8/C9 family, present in serum. C8 is a terminal component of the complement system, part of both the complement membrane attack complex (MAC), and important to MAC assembly. C8 binds to the C5b-7 complex, anchored to the membrane, creating C5b-8. C5b-8 binds C9 and catalyses the polymerization of C9 molecules to form C5b-9 (MAC). C8 is thought to contain lipid binding sites, facilitating the insertion of MAC into the membrane. Defects in the alpha chain of C8 can result in complement C8 deficiency type I. Furthermore, C8 deficiencies can cause recurring bacterial infections, in particular from <i>Neisseria meningitides</i>.
Histology Positive Control Tissue	Kidney from post streptococcal glomerulonephritis patients
References	1. Rosa, D.D. <i>et al.</i> (2004) Deficiency of the eighth component of complement associated with recurrent meningococcal meningitis--case report and literature review. Braz J Infect Dis. 8 (4): 328-30. 2. Schreck, S.F. <i>et al.</i> (2000) Human complement protein C8 gamma. Biochim Biophys Acta. 1482 (1-2): 199-208.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10334 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10334.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP

Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Human Anti Mouse IgG1 (HCA036...)	HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

Recommended Useful Reagents

[HISTAR DETECTION SYSTEM \(STAR3000A\)](#)

[HISTAR DETECTION SYSTEM \(STAR3000B\)](#)

[HISTAR DETECTION SYSTEM \(STAR3000C\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M311681:171215'

Printed on 16 Dec 2017