

## Datasheet: MCA340PE

<b>Description:</b>	MOUSE ANTI RAT CD45RA:RPE
<b>Specificity:</b>	CD45RA (B CELLS ONLY)
<b>Format:</b>	RPE
<b>Product Type:</b>	Monoclonal Antibody
<b>Clone:</b>	OX-33
<b>Isotype:</b>	IgG1
<b>Quantity:</b>	100 TESTS

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			Neat - 1/10

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rat		
<b>Product Form</b>	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
<b>Reconstitution</b>	Reconstitute with 1 ml distilled water		
<b>Max Ex/Em</b>	<b>Fluorophore</b>	<b>Excitation Max (nm)</b>	<b>Emission Max (nm)</b>
	RPE 488nm laser	496	578
<b>Buffer Solution</b>	Phosphate buffered saline		
<b>Preservative</b>	0.09% Sodium Azide		
<b>Stabilisers</b>	1% Bovine Serum Albumin		
	5% Sucrose		
<b>Immunogen</b>	Purified Rat spleen L-CA		
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P04157</a> <a href="#">Related reagents</a>		
	<b>Entrez Gene:</b> <a href="#">24699</a> Ptprc <a href="#">Related reagents</a>		

<b>Fusion Partners</b>	Spleen cells from immunised BALB/c mice were fused with cells of the NSO/U mouse myeloma cell line.
<b>Specificity</b>	<b>Mouse anti Rat CD45RA antibody, clone OX-33</b> is directed against a high molecular weight band of the leucocyte common antigen. MRC OX-33 only labels B-cells among thoracic duct lymphocytes with little labelling in bone marrow and none on thymocytes ( <a href="#">Barclay et al. 1987</a> ).
<b>Flow Cytometry</b>	Use 10ul of the suggested working dilution to label 10 <sup>6</sup> cells in 100ul.
<b>References</b>	<ol style="list-style-type: none"> <li>1. Woollett, G.R. <i>et al.</i> (1985) Molecular and antigenic heterogeneity of the rat leukocyte-common antigen from thymocytes and T and B lymphocytes. <a href="#">Eur J Immunol. 15 (2): 168-73.</a></li> <li>2. Barclay, A.N. <i>et al.</i> (1987) Lymphocyte specific heterogeneity in the rat leucocyte common antigen (T200) is due to differences in polypeptide sequences near the NH2-terminus. <a href="#">EMBO J. 6:1259-1264.</a></li> <li>3. Barclay, A.N. <i>et al.</i> (1981) The localization of populations of lymphocytes defined by monoclonal antibodies in rat lymphoid tissues. <a href="#">Immunology. 42:593-600.</a></li> <li>4. Whiteland, J.L. <i>et al.</i> (1995) Immunohistochemical detection of T-cell subsets and other leukocytes in paraffin-embedded rat and mouse tissues with monoclonal antibodies. <a href="#">J Histochem Cytochem. 43 (3): 313-20.</a></li> <li>5. McCall, M.N. <i>et al.</i> (1992) Expression of soluble isoforms of rat CD45. Analysis by electron microscopy and use in epitope mapping of anti-CD45R monoclonal antibodies. <a href="#">Immunology. 76: 310-7.</a></li> <li>6. Cho, W.S. <i>et al.</i> (2012) NiO and Co3O4 nanoparticles induce lung DTH-like responses and alveolar lipoproteinosis. <a href="#">Eur Respir J. 39 (3): 546-57.</a></li> <li>7. Williamson Lauren L. <i>et al.</i> (2015) Got worms? Perinatal exposure to helminths prevents persistent immune sensitization and cognitive dysfunction induced by early-life infection <a href="#">Brain, Behavior, and Immunity. pii: S0889-1591(15)00240-8.</a></li> <li>8. Pilatz, A. <i>et al.</i> (2015) Experimental <i>Escherichia coli</i> epididymitis in rats: assessment of testicular involvement in a long-term follow-up. <a href="#">Andrologia. 47 (2): 160-7.</a></li> <li>9. Stefanski, V. <i>et al.</i> (2013) Differential effect of severe and moderate social stress on blood immune and endocrine measures and susceptibility to collagen type II arthritis in male rats. <a href="#">Brain Behav Immun. 29: 156-65.</a></li> <li>10. Trama, A.M. <i>et al.</i> (2012) Lymphocyte phenotypes in wild-caught rats suggest potential mechanisms underlying increased immune sensitivity in post-industrial environments. <a href="#">Cell Mol Immunol. 9 (2): 163-74.</a></li> <li>11. Okamura, T. <i>et al.</i> (2013) Phenotypic Characterization of LEA Rat: A New Rat Model of Nonobese Type 2 Diabetes. <a href="#">J Diabetes Res. 2013: 986462.</a></li> <li>12. Denecke, C. <i>et al.</i> (2013) Synergistic effects of prolonged warm ischemia and donor age on the immune response following donation after cardiac death kidney transplantation. <a href="#">Surgery. 153 (2): 249-61.</a></li> <li>13. Marolda, R. <i>et al.</i> (2013) Differential targeting of immune-cells by Pixantrone in experimental myasthenia gravis. <a href="#">J Neuroimmunol. 258 (1-2): 41-50.</a></li> <li>14. Han, X. <i>et al.</i> (2013) <i>Porphyromonas gingivalis</i> infection-associated periodontal bone resorption is dependent on receptor activator of NF-κB ligand. <a href="#">Infect Immun. 81 (5): 1502-9.</a></li> <li>15. Pongratz, G. <i>et al.</i> (2015) A sustained high fat diet for two years decreases IgM and IL-1 beta in ageing Wistar rats. <a href="#">Immun Ageing. 12: 12.</a></li> <li>16. Lu, J.H. <i>et al.</i> (2015) GABAergic neurons in cerebellar interposed nucleus modulate cellular and humoral immunity via hypothalamic and sympathetic pathways. <a href="#">J Neuroimmunol. 283: 30-8.</a></li> </ol>
<b>Storage</b>	<p>Prior to reconstitution store at +4°C. Following reconstitution store at +4°C.</p> <p>DO NOT FREEZE.</p>

This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

---

<b>Shelf Life</b>	18 months from date of reconstitution
-------------------	---------------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10075 available at: 10075: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf</a>
--------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

## Related Products

### Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA1209PE\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

'M300903:170106'

**Printed on 17 May 2018**

---

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)