

Datasheet: MCA424B

Description:	RAT ANTI MOUSE IgG HEAVY CHAIN:Biotin
Specificity:	IgG HEAVY CHAIN
Format:	Biotin
Product Type:	Monoclonal Antibody
Clone:	COCKTAIL
Isotype:	Cocktail
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			1.0ug/ml
Western Blotting			■	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG conjugated to Biotin - liquid
Preparation	Purified IgG prepared by affinity chromatography on mouse Ig/Sepharose 4B from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative	0.1% Sodium Azide
Stabilisers	50% Glycerol
Approx. Protein Concentrations	IgG concentration 1 mg/ml
Immunogen	Purified IgG from BALB/c mice.
External Database Links	UniProt: P01869 Related reagents

[P01865](#) [Related reagents](#)
[P03987](#) [Related reagents](#)
[P01867](#) [Related reagents](#)
[P01864](#) [Related reagents](#)
[P01868](#) [Related reagents](#)

Entrez Gene:

[16017](#) Ighg1 [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)
[16016](#) Ighg2b [Related reagents](#)
[16017](#) Ighg1 [Related reagents](#)
[380793](#) Igh-1a [Related reagents](#)
[380795](#) AI324046 [Related reagents](#)

Synonyms	Igh-4
Fusion Partners	Spleen cells from immunised LOU/c rats were fused with cells of the IR983F rat myeloma cell line.
Specificity	Rat anti Mouse IgG Heavy Chain antibody cocktail has been shown to be specific for the gamma heavy chain of mouse immunoglobulin by immunoblotting techniques. It recognizes IgG2a and 2b strongly and IgG1 and IgG3 more weakly.
References	1. Bazin, H. (1982) Production of rat monoclonal antibodies with the LOU rat non secreting IR983F myeloma cell line. Prot. Biol. Fluids. Peters Ed: 29th Colloquim (1981). Pergamon Press, Oxford and N.Y., pp 615-618. 2. Bazin, H. <i>et al.</i> (1984) Rat monoclonal antibodies. II. A rapid and efficient method of purification from ascitic fluid or serum J. Immunol. Methods 71: 9 -16
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10328 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10328.pdf
Regulatory	For research purposes only

North & South America Tel: +1 800 265 7376
 Fax: +1 919 878 3751
 Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
 Fax: +44 (0)1865 852 739
 Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
 Fax: +49 (0) 89 8090 95 50
 Email: antibody_sales_de@bio-rad.com

'M314221:180412'

Printed on 13 Apr 2018

