

Datasheet: MCA4699GA

Description:	MOUSE ANTI HUMAN CD266
Specificity:	CD266
Other names:	TWEAK-R
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	ITEM-4
Isotype:	IgG2b
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/25 - 1/200
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA			■	
Immunoprecipitation			■	
Western Blotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0mg/ml

Immunogen	Human CD266 transfected P815 cells.
External Database Links	UniProt: Q9NP84 Related reagents Entrez Gene: 51330 TNFRSF12A Related reagents
Synonyms	FN14
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the P3U1 mouse myeloma cell line.
Specificity	Mouse anti Human CD266 antibody, clone ITEM-4 recognizes human Tumor necrosis factor receptor superfamily member 12A, also known as CD266, Fibroblast growth factor-inducible immediate-early response protein 14 or Tweak-receptor. CD266 is a 102 amino acid ~14 kDa type I transmembrane glycoprotein and member of the tumour necrosis factor receptor superfamily (TNFRSF12A) expressed by a variety of tissues, including the heart, lung and kidney, and at elevated levels by human brain and liver tumours. The interaction between CD266 and its ligand CD255 (TWEAK), induces the proliferation and migration of endothelial cells and also angiogenesis, and studies have shown that CD255-induced apoptosis of some tumour cell lines is mediated by CD266 (Nakayama et al. 2003). Mouse anti Human CD266 antibody, clone ITEM-4 is cross-reactive with mouse CD266, and also with Rat-1 fibroblastoid cells in flow cytometry.
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
References	1. Nakayama, M. <i>et al.</i> (2003) Fibroblast growth factor-inducible 14 mediates multiple pathways of TWEAK-induced cell death. J Immunol. 170 (1): 341-8. 2. Nakayama, M. <i>et al.</i> (2003) Characterization of murine TWEAK and its receptor (Fn14) by monoclonal antibodies. Biochem Biophys Res Commun. 306 (4): 819-25. 3. Winkles, J.A. and Brown, S. (2010) Targeted Therapy of Fn14-Positive Breast Tumors Using a TWEAK-Cytotoxin Fusion Protein or Noncovalent Complex. U.S. Army Medical Research and Materiel Command, Annual Report 4. Girgenrath, M. <i>et al.</i> (2006) TWEAK, via its receptor Fn14, is a novel regulator of mesenchymal progenitor cells and skeletal muscle regeneration. EMBO J. 25: 5826-39.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG2b (HCA038...)	FITC , HRP

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL \(MCA691\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M316007:180503'

Printed on 05 May 2018