

Datasheet: MCA4758

Description:	MOUSE ANTI PSEUDOMONAS AERUGINOSA SEROTYPE 8
Specificity:	PSEUDOMONAS AERUGINOSA SEROTYPE 8
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	1017/87
Isotype:	IgM
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			
Immunoassay	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Ig fraction - liquid
Preparation	Ig fraction prepared by ammonium sulphate precipitation from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgM concentration 1.0 mg/ml
Immunogen	<i>P. aeruginosa</i> serotype 8 (ATCC 33355).

Specificity **Mouse anti *Pseudomonas aeruginosa* serotype 8 antibody, clone 1017/87** recognizes *P. aeruginosa* serotype 8. as described in ATCC:33355.

P. aeruginosa is a gram-negative, opportunistic pathogen found in soil and water. It is a major cause of hospital acquired infection, and can affect any part of the body including the urinary, respiratory and gastrointestinal tracts in addition to skin, soft tissue, bones and joints.

Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgM (STAR86...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgM (STAR138...)	Alk. Phos.
Human Anti Mouse IgM (HCA040...)	FITC , HRP

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	--	------------------	--	---------------	--

'M317186:180619'

Printed on 20 Jun 2018