

Datasheet: MCA4755EL

Description:	MOUSE ANTI HUMAN CD205:Low Endotoxin
Specificity:	CD205
Other names:	DEC205
Format:	Low Endotoxin
Product Type:	Monoclonal Antibody
Clone:	HD30
Isotype:	IgG1
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/10 - 1/200
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA			■	
Immunoprecipitation			■	
Western Blotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Species Cross Reactivity	Reacts with: Monkey Does not react with: Mouse N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	<0.01EU/ug

Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Purified extracellular domain of human CD205-IgG1 Fc fusion protein.
External Database Links	UniProt: O60449 Related reagents Entrez Gene: 4065 LY75 Related reagents
Synonyms	CD205, CLEC13B
Fusion Partners	Spleen cells from immunised mice were fused with cells of the mouse SP2/0-Ag14 myeloma cell line.
Specificity	Mouse anti Human CD205 antibody, clone HD30 recognizes human Lymphocyte antigen 75, otherwise known as CD205, DEC205, C-type lectin domain family 13 member B or gp200-MR6. CD205 is a 1695 amino acid ~205 kDa single pass type I transmembrane glycoprotein and member of the macrophage mannose receptor family of C-type lectins, predominantly expressed by human dendritic cells. CD205 is a recycled endocytic receptor involved in directing antigen uptake for processing and presentation by MHC class I. Multiple isoforms of human CD205, generated by alternative splicing have been described. Mouse anti Human CD205 antibody, clone HD30 recognizes the canonical isoform 4, reactivity with other isoforms has not been evaluated but both isoforms 2 and 3 contain most of the immunogen sequence used to generate clone HD30 making reactivity likely. However the antibody is not expected to be reactive with isoforms 1 and 5 which lack the immunogen sequence (UniProt : 060449).
Flow Cytometry	Use 10ul of suggested working dilution 10 ⁶ cells in 100ul.
Histology Positive Control Tissue	Human or monkey spleen and lymph node.
References	1. Park, C.G. <i>et al.</i> (2012) Generation of anti-human DEC205/CD205 monoclonal antibodies that recognize epitopes conserved in different mammals. J Immunol Methods. 377 (1-2): 15-22.
Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10162 available at: 10162: https://www.bio-rad-antibodies.com/uploads/MSDS/10162.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG1 (HCA036...)	HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:Low Endotoxin \(MCA928EL\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M316023:180503'

Printed on 05 May 2018