

Datasheet: MCA4755PE

Description:	MOUSE ANTI HUMAN CD205:RPE
Specificity:	CD205
Other names:	DEC205
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	HD30
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human								
Species Cross Reactivity	Reacts with: Monkey Does not react with: Mouse N.B. Antibody reactivity and working conditions may vary between species.								
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized								
Reconstitution	Reconstitute with 1ml distilled water								
Max Ex/Em	<table><tr><th>Fluorophore</th><th>Excitation Max (nm)</th><th>Emission Max (nm)</th></tr><tr><td>RPE 488nm laser</td><td>496</td><td>578</td></tr></table>	Fluorophore	Excitation Max (nm)	Emission Max (nm)	RPE 488nm laser	496	578		
Fluorophore	Excitation Max (nm)	Emission Max (nm)							
RPE 488nm laser	496	578							
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant								
Buffer Solution	Phosphate buffered saline								
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin								
Immunogen	Purified extracellular domain of human CD205-IgG1 Fc fusion protein.								

**External Database
Links**

UniProt:

[O60449](#)

[Related reagents](#)

Entrez Gene:

[4065](#)

LY75

[Related reagents](#)

Synonyms

CD205, CLEC13B

Fusion Partners

Spleen cells from immunised mice were fused with cells of the mouse SP2/0-Ag14 myeloma cell line.

Specificity

Mouse anti Human CD205 antibody, clone HD30 recognizes human Lymphocyte antigen 75, otherwise known as CD205, DEC205, C-type lectin domain family 13 member B or gp200-MR6. CD205 is a 1695 amino acid ~205 kDa single pass type I transmembrane glycoprotein and member of the macrophage mannose receptor family of C-type lectins, predominantly expressed by human dendritic cells.

CD205 is a recycled endocytic receptor involved in directing antigen uptake for processing and presentation by MHC class I. Multiple isoforms of human CD205, generated by alternative splicing have been described. Mouse anti Human CD205 antibody, clone HD30 recognizes the canonical isoform 4, reactivity with other isoforms has not been evaluated but both isoforms 2 and 3 contain most of the immunogen sequence used to generate clone HD30 making reactivity likely. However the antibody is not expected to be reactive with isoforms 1 and 5 which lack the immunogen sequence ([UniProt : 060449](#)).

Flow Cytometry

Use 10ul of suggested working dilution 10⁶ cells in 100ul.

References

1. Park, C.G. *et al.* (2012) Generation of anti-human DEC205/CD205 monoclonal antibodies that recognize epitopes conserved in different mammals. [J Immunol Methods. 377 \(1-2\): 15-22.](#)

Storage

Prior to reconstitution store at +4°C.

After reconstitution store at +4°C.

DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

12 months from date of reconstitution.

**Health And Safety
Information**

Material Safety Datasheet documentation #10075 available at:
10075: <https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf>

Regulatory

For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA928PE\)](#)

Recommended Useful Reagents

[HUMAN SEROBLOCK \(BUF070A\)](#)

[HUMAN SEROBLOCK \(BUF070B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com 'M300727:170106'	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

Printed on 05 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)