

Datasheet: MCA4771

Description:	MOUSE ANTI HUMAN TALIN-2
Specificity:	TALIN-2
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	53.8
Isotype:	IgG2b
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin	■			1/25 - 1/100
ELISA			■	
Immunoprecipitation	■			
Western Blotting	■			
Immunofluorescence (1)	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1) **Methanol treatment is required prior to antibody staining**

Target Species	Human
Species Cross Reactivity	Reacts with: Mouse, Rabbit N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml

Immunogen	Recombinant Mouse talin-2 (aa2300-2542).
External Database Links	UniProt: Q9Y4G6 Related reagents Entrez Gene: 83660 TLN2 Related reagents
Synonyms	KIAA0320
Specificity	Mouse anti Human Talin-2 antibody, clone 53.8 (also designated 121A) recognizes Talin-2, a high molecular weight skeletal protein involved in focal adhesion plaques. It plays a role in linking integrin to the actin cytoskeleton. It interacts with phosphatidylinositol-4-phosphate 5-kinase type I gamma (PIP5K1C), inducing its kinase activity. Mouse anti Human Talin-2 antibody, clone 53.8 recognizes an epitope within the range of amino acids 2470-2500.
Histology Positive Control Tissue	Tonsil
References	1. Liu, J. <i>et al.</i> (2015) Talin determines the nanoscale architecture of focal adhesions. Proc Natl Acad Sci U S A. 112 (35): E4864-73. 2. Prince, H.E. & Yeh, C. (2013) Reactivity of human IgM binding murine monoclonal 6B6C1 (IgG2a) with other murine monoclonal IgG antibodies. J Clin Lab Anal. 27 (1): 27-30. 3. Qi, L. <i>et al.</i> (2016) Talin2-mediated traction force drives matrix degradation and cell invasion. J Cell Sci. 129 (19): 3661-74. 4. Ichikawa, Y. <i>et al.</i> (2017) Modulation of caveolins, integrins and plasma membrane repair proteins in anthracycline-induced heart failure in rabbits. PLoS One. 12 (5): e0177660.
Further Reading	1. Critchley, D.R. & Gingras, A.R. (2008) Talin at a glance. J Cell Sci. 121 (Pt 9): 1345-7. 2. Monkley SJ <i>et al.</i> (2001) Analysis of the mammalian talin2 gene TLN2. Biochem Biophys Res Commun. 286 (5): 880-5.
Storage	Store at +4°C or at -20°C if preferred. This product should be stored undiluted. Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG2b (HCA038...)	FITC , HRP

Recommended Negative Controls

[MOUSE IgG2b NEGATIVE CONTROL \(MCA691\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M316030:180503'

Printed on 05 May 2018