

Datasheet: MCA494GA

Description:	MOUSE ANTI RAT CD25
Specificity:	CD25
Other names:	IL-2R ALPHA CHAIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	NDS 61
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/50 - 1/100
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA	■			
Immunoprecipitation			■	
Western Blotting			■	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1 mg/ml
Immunogen	Rat interleukin-2 receptor.

**External Database
Links**

UniProt:

[P26897](#)

[Related reagents](#)

Entrez Gene:

[25704](#)

IL2ra

[Related reagents](#)

Fusion Partners

Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS1 myeloma cell line.

Specificity

Mouse anti Rat CD25 antibody, clone NDS 61 recognizes the rat CD25 cell surface antigen, a ~55 kDa molecule that is the low affinity alpha chain of the IL-2 receptor. Mouse anti Rat CD25 antibody, clone NDS 61 recognizes a different epitope on CD25 to Mouse anti Rat CD25, clone OX-39 ([MCA273R](#)) .

Mouse anti Rat CD25 antibody, clone NDS 61 is reported to block IL-2 induced cell proliferation ([Tellides *et al.* 1989](#)).

Flow Cytometry

Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

1. Tellides, G. *et al.* (1987) Functional blocking of the interleukin-2 receptor (IL-2R) may be important in the efficacy of IL-2R antibody therapy. [Transplant Proc. 19 \(5\): 4231-3.](#)
2. Tellides, G. *et al.* (1989) Mechanism of action of interleukin-2 receptor (IL-2R) monoclonal antibody (MAb) therapy: target cell depletion or inhibition of function? [Transplant Proc. 21 \(1 Pt 1\): 997-8.](#)
3. Dallman, M.J. *et al.* (1988) Cytotoxicity and IL 2 reactivity of cells from rejected or nonrejected allografts. [Transplant Proc. 20 \(2\): 226-8.](#)
4. Bjersing, J.L. *et al.* (2002) Loss of ileal IgA+ plasma cells and of CD4+ lymphocytes in ileal Peyer's patches of vitamin A deficient rats. [Clin Exp Immunol. 130: 404-8.](#)
5. Gassel, H.J. *et al.* (2000) Tolerance of rat liver allografts induced by short-term selective immunosuppression combining monoclonal antibodies directed against CD25 and CD54 with subtherapeutic cyclosporine. [Transplantation. 69: 1058-67.](#)
6. Herías, M.V. *et al.* (1999) Immunomodulatory effects of *Lactobacillus plantarum* colonizing the intestine of gnotobiotic rats. [Clin Exp Immunol. 116: 283-90.](#)
7. Penttilä, I.A. *et al.* (2003) Maternal milk regulation of cell infiltration and interleukin 18 in the intestine of suckling rat pups. [Gut. 52: 1579-86.](#)
8. Sawitzki, B. *et al.* (2004) IFN-gamma regulation in anti-CD4 antibody-induced T cell unresponsiveness. [J Am Soc Nephrol. 15: 695-703.](#)
9. Wu, H. *et al.* (2005) DNA vaccination with naked DNA encoding MCP-1 and RANTES protects against renal injury in adriamycin nephropathy. [Kidney Int. 67: 2178-86.](#)
10. Stanisavljević S *et al.* (2016) Gut-associated lymphoid tissue, gut microbes and susceptibility to experimental autoimmune encephalomyelitis. [Benef Microbes. Feb 3 \[Epub ahead of print\]](#)

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

Health And Safety

Material Safety Datasheet documentation #10040 available at:

Information 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG1 (HCA036...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA1209\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M319170:180719'

Printed on 01 Aug 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)