

Datasheet: MCA555

Description:	MOUSE ANTI HUMAN CD43
Specificity:	CD43
Other names:	LEUKOSIALIN
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	DFT-1
Isotype:	IgG1
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/100
Immunohistology - Frozen	■			
Immunohistology - Paraffin	■			1/10 - 1/100
ELISA			■	
Immunoprecipitation			■	
Western Blotting			■	
Immunofluorescence	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1 mg/ml
Immunogen	KG-1 cells.

**External Database
Links**

UniProt:

[P16150](#) [Related reagents](#)

Entrez Gene:

[6693](#) SPN [Related reagents](#)

Synonyms

CD43

Specificity

Mouse anti Human CD43 antibody, clone DFT-1 recognizes human CD43 also known as leukosialin or sialophorin. CD43 is a 381 amino acid single pass type I transmembrane glycoprotein present on the cell surface of multiple cell types including thymocytes, T-lymphocytes, neutrophils, plasma cells and platelets. It is lacking on the surface of most B-lymphocytes. CD34 is a major membrane sialoglycoprotein expressed on T lymphocytes as two predominant glycoforms, a 115kDa form present on all T-cells and a 130kDa form indicative of T-cell activation ([Barren *et al.* 1997](#)) and acts as a T-cell counter-receptor for [CD169/Siglec-1](#) ([van den Berg *et al.* 2001](#)). Clone DFT-1 binds to both glycoforms of human CD43 ([Kyoizumi *et al.* 2004](#))

CD34 expression is deficient and/or defective in patients suffering from the X-linked chromosomal immunodeficiency condition [Wiskott-Aldrich syndrome](#) and acts as a ligand for [CD54](#) ([Rosenstein *et al.* 1991](#)).

Mouse anti human CD43, clone DFT-1 has been used successfully for the evaluation of CD43 expression in patients with Myelodysplastic syndromes ([Kyriakou *et al.* 2006](#))

Flow Cytometry

Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

Immunohistology

This product does not require protein digestion pre-treatment of paraffin sections however, staining of paraffin embedded formalin fixed tissues may be improved in some cases by enzyme pretreatment using Trypsin.

**Histology Positive
Control Tissue**

Tonsil

References

1. Stross, W.P. *et al.* (1989) Molecule detected in formalin fixed tissue by antibodies MT1, DF-T1, and L60 (Leu-22) corresponds to CD43 antigen. [J Clin Pathol. 42 \(9\): 953-61.](#)
2. Kojima N *et al.* (1994) Role of cell surface O-linked oligosaccharides in adhesion of HL60 cells to fibronectin: regulation of integrin-dependent cell adhesion by O-linked oligosaccharide elongation. [Exp Cell Res. 214 \(2\): 537-42.](#)
3. Anzai, N. *et al.* (1999) Modulation of integrin function in hematopoietic progenitor cells by CD43 engagement: possible involvement of protein tyrosine kinase and phospholipase C-gamma. [Blood. 93 \(10\): 3317-26.](#)
4. Matsuo, A. *et al.* (1996) Expression of CD43 in human microglia and its downregulation in Alzheimer's disease. [J Neuroimmunol. 71 \(1-2\): 81-6.](#)
5. Nonomura, C. *et al.* (2008) CD43, but not P-selectin glycoprotein ligand-1, functions as an E-selectin counter-receptor in human pre-B-cell leukemia NALL-1. [Cancer Res. 68 \(3\): 790-9.](#)

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend

microcentrifugation before use.

Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG1 (HCA036...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M319191:180719'

Printed on 01 Aug 2018