

Datasheet: MCA555PE

Description:	MOUSE ANTI HUMAN CD43:RPE
Specificity:	CD43
Other names:	LEUKOSIALIN
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	DFT-1
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Human		
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
Reconstitution	Reconstitute with 1 ml distilled water		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578
Preparation	Purified IgG prepared by affinity chromatography on Protein G from		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1%	Bovine Serum Albumin	
	5%	Sucrose	
Immunogen	KG-1 cells.		
External Database Links	UniProt: P16150 Related reagents		

Entrez Gene:[6693](#) SPN [Related reagents](#)

Synonyms	CD43
Specificity	Mouse anti Human CD43 antibody, clone DFT-1 recognizes human CD43 also known as leukosialin or sialophorin. CD43 is a 381 amino acid single pass type I transmembrane glycoprotein present on the cell surface of multiple cell types including thymocytes, T-lymphocytes, neutrophils, plasma cells and platelets. It is lacking on the surface of most B-lymphocytes. CD34 is a major membrane sialoglycoprotein expressed on T lymphocytes as two predominant glycoforms, a 115kDa form present on all T-cells and a 130kDa form indicative of T-cell activation (Barren et al. 1997) and acts as a T-cell counter-receptor for CD169/Siglec-1 (van den Berg et al. 2001). Clone DFT-1 binds to both glycoforms of human CD43 (Kyoizumi et al. 2004) CD34 expression is deficient and/or defective in patients suffering from the X-linked chromosomal immunodeficiency condition Wiskott-Aldrich syndrome and acts as a ligand for CD54 (Rosenstein et al. 1991). Mouse anti human CD43, clone DFT-1 has been used successfully for the evaluation of CD43 expression in patients with Myelodysplastic syndromes (Kyriakou et al. 2006)
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells or 100ul whole blood
References	1. Stross, W.P. <i>et al.</i> (1989) Molecule detected in formalin fixed tissue by antibodies MT1, DF-T1, and L60 (Leu-22) corresponds to CD43 antigen. J Clin Pathol. 42 (9): 953-61. 2. Kojima N <i>et al.</i> (1994) Role of cell surface O-linked oligosaccharides in adhesion of HL60 cells to fibronectin: regulation of integrin-dependent cell adhesion by O-linked oligosaccharide elongation. Exp Cell Res. 214 (2): 537-42. 3. Anzai, N. <i>et al.</i> (1999) Modulation of integrin function in hematopoietic progenitor cells by CD43 engagement: possible involvement of protein tyrosine kinase and phospholipase C-gamma. Blood. 93 (10): 3317-26. 4. Matsuo, A. <i>et al.</i> (1996) Expression of CD43 in human microglia and its downregulation in Alzheimer's disease. J Neuroimmunol. 71 (1-2): 81-6. 5. Nonomura, C. <i>et al.</i> (2008) CD43, but not P-selectin glycoprotein ligand-1, functions as an E-selectin counter-receptor in human pre-B-cell leukemia NALL-1. Cancer Res. 68 (3): 790-9.
Storage	Prior to reconstitution store at +4°C. Following reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet documentation #10075 available at: 10075: https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

MOUSE IgG1 NEGATIVE CONTROL:RPE (MCA928PE)

Recommended Useful Reagents

HUMAN SEROBLOCK (BUF070A)

HUMAN SEROBLOCK (BUF070B)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M300883:170106'

Printed on 01 Aug 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)