

Datasheet: MCA5799GA

Description:	MOUSE ANTI RAT CD1d
Specificity:	CD1d
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	WTH1
Isotype:	IgG2a
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/2,000
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA	■			
Immunoprecipitation	■			
Western Blotting (1)	■			
Functional Assays (2)	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

(1)**Clone WTH1 recognizes CD1d under non-reducing conditions, see [Monzon-Casanova, E. et al.](#) for details.**

(2)This antibody contains sodium azide, removal by dialysis is recommended prior to use in functional assays. Bio-Rad recommend the use of EQU003 for this purpose.

Target Species	Rat
Species Cross Reactivity	Reacts with: Mouse Does not react with: Cotton rat N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.02% Sodium Azide (NaN ₃)

Approx. Protein Concentrations	IgG concentration 0.5mg/ml
Immunogen	Rat CD1d-transduced M12.4.1.C3 cells.
External Database Links	UniProt: P11609 Related reagents Q63493 Related reagents Entrez Gene: 12479 Cd1d1 Related reagents 25109 Cd1d1 Related reagents
Synonyms	Cd1.1, Cd1d1
Specificity	Mouse anti Rat CD1d antibody, clone WTH1 recognizes rat CD1d, a glycoprotein expressed on cortical thymocytes, certain T-cell leukaemias and various other antigen presenting cells. CD1d binds self and non-self glycolipids and presents them to T-cell receptors on natural killer cells. When activated, natural killer cells produce Th1 and Th2 cytokines. Clone WTH1 recognizes a different, non-overlapping epitope to clone WTH2 and binds with high affinity to the mature CD1d heavy chain independently of beta-2-microglobulin. Mouse anti Rat CD1d antibody, clone WTH1 efficiently blocks rat and mouse invariant natural killer T cell (iNKT cell) response to a-GalCer (an antigen for iNKT cells) and has a slightly higher affinity to native CD1d than clone WTH2. Clone WTH1 cross reacts weakly with mouse tissue in IHC and western blotting, but reacts strongly with both rat and mouse antigens in flow cytometry and immunoprecipitation.
Flow Cytometry	Use 10ul of the suggested working dilution to label 1x10 ⁶ cells in 100ul.
References	1. Monzon-Casanova, E. <i>et al.</i> (2010) CD1d expression in paneth cells and rat exocrine pancreas revealed by novel monoclonal antibodies which differentially affect NKT cell activation. PLoS One 5:e13089 2. Fichtner AS <i>et al.</i> (2015) Function and expression of CD1d and invariant natural killer T-cell receptor in the cotton rat (<i>Sigmodon hispidus</i>). Immunology. 146 (4): 618-29. 3. Monzon-Casanova, E. <i>et al.</i> (2013) Direct identification of rat iNKT cells reveals remarkable similarities to human iNKT cells and a profound deficiency in LEW rats. Eur J Immunol. 43 (2): 404-15.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Human Anti Mouse IgG2a (HCA037...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR13...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG2a NEGATIVE CONTROL \(MCA1210\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M319204:180719'

Printed on 01 Aug 2018