

Datasheet: MCA5936GA

Description:	MOUSE ANTI CHICKEN CD34
Specificity:	CD34
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	AV138
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA	■			
Immunoprecipitation			■	
Western Blotting	■			
Immunocytochemistry (1)	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

(1) It has been demonstrated that this antibody can be used to detect chicken CD34 expressed in COS7 cells using an AEC substrate detection system.

Target Species	Chicken
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0mg/ml
Immunogen	Chicken CD34/hulg fusion protein.

**External Database
Links**

UniProt:

[E1BUT3](#)

[Related reagents](#)

Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the Sp2/0 Ag14 Myeloma cell line.
Specificity	Mouse anti Chicken CD34 antibody, clone AV138 recognizes chicken CD34, a cell surface antigen expressed by, endothelial progenitor cells, endothelial cells and haematopoietic stem cells. Mouse anti Chicken CD34 antibody, clone AV138 is one of a range of antibodies that have been developed by the Roslin Institute as part of the Avian immunology toolbox project. In addition to this anti chicken CD34 antibody, other antibodies to chicken CD antigens and cytokines from this project are now available from Bio-Rad.
ELISA	This product is suitable for use in indirect ELISA applications.
Western Blotting	Clone AV138 detects a band of approximately 70 - 80kDa on COS cell supernatant transfected with chicken CD34 (extracellular)/V5HIS6 tagged protein under reducing conditions.
Further Reading	1. Cairns, D.M. <i>et al.</i> (2012) Interplay of Nkx3.2, Sox9 and Pax3 regulates chondrogenic differentiation of muscle progenitor cells. PLoS One. 7: e39642.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch.
Acknowledgements	This reagent was produced by the WellcomeTrust-funded National Avian Resource Facility (Grant number WT099164MA).
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP

Human Anti Mouse IgG1 (HCA036...) [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®549](#),
[DyLight®649](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M317297:180625'

Printed on 25 Jun 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)