

Datasheet: MCA5970GA

Description:	MOUSE ANTI HUMAN OSTEOPROTEGERIN
Specificity:	OSTEOPROTEGERIN
Other names:	TNFRSF11B
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	E01-10F11
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			
Western Blotting	■			1/1000

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from ascites
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	<i>E. coli</i> derived recombinant protein corresponding to aa 22-201 of human osteoprotegerin
External Database Links	UniProt: O00300 Related reagents Entrez Gene: 4982 TNFRSF11B Related reagents

Synonyms	OCIF, OPG
Fusion Partners	Spleen cells from immunised BALB/c mice were fused with cells of the mouse SP2/0 myeloma cell line
Specificity	Mouse anti Human osteoprotegerin, clone E01-10F11, recognizes osteoprotegerin (OPG), a cytokine receptor and member of the tumour necrosis factor (TNF) receptor superfamily. The role of OPG involves it acting as a decoy receptor for the ligand of CD254 and thereby preventing its function in osteoclastogenesis. In addition to this it may act as a decoy receptor for TNF related apoptosis inducing ligand (TRAIL) (Emery et al. 1998 and Schoppet et al. 2002). The combination of these two roles allows OPG to inhibit NFkB's regulatory effects on inflammation, skeletal, and vascular systems. OPG is highly expressed in a wide range of tissues, including adult heart, lung, liver, prostate, ovary and bone marrow. More recently a role for OPG in cardiovascular disease has emerged and it is now clinically relevant biomarker of vascular health (Bierre, 2013).
Further Reading	1. Emery, J.G. <i>et al.</i> (1998) Osteoprotegerin is a receptor for the cytotoxic ligand TRAIL. J Biol Chem. 273 (23): 14363-7. 2. Schoppet, M. <i>et al.</i> (2004) Localization of osteoprotegerin, tumor necrosis factor-related apoptosis-inducing ligand, and receptor activator of nuclear factor-kappaB ligand in Mönckeberg's sclerosis and atherosclerosis. J Clin Endocrinol Metab. 89 (8): 4104-12. 3. Bjerre, M. (2013) Osteoprotegerin (OPG) as a biomarker for diabetic cardiovascular complications. Springerplus. 2: 658.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800

Goat Anti Mouse IgG (STAR70...) [FITC](#)

Rabbit Anti Mouse IgG (STAR13...) [HRP](#)

Human Anti Mouse IgG1 (HCA036...) [HRP](#)

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA928\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M261641:140925'

Printed on 01 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)