

Datasheet: MCA6011GA

Description:	MOUSE ANTI CHICKEN CD268 / BAFF-R
Specificity:	CD268 / BAFF-R
Other names:	TNFRSF 13C
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	2C4
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			
Immunohistology - Frozen	■			
Immunohistology - Paraffin	■			
ELISA			■	
Immunoprecipitation			■	
Western Blotting			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Chicken
Species Cross Reactivity	Reacts with: Turkey, Quail, Duck N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Chicken BAFF-R transfected HEK293 cells

Specificity	Mouse anti Chicken CD268 / BAFF-R, clone 2C4, recognises the chicken homologue of human BAFF-R, known also as CD268 and TNFRSF 13C, a cell surface antigen expressed by avian B-cells. Chicken BAFF-R (B cell activating factor receptor) differs from mammalian BAFF-R in that it is produced in B-cells (Reddy et al. 2008) and is not expressed by chicken T-cells. In addition to chicken, clone 2C4 has been shown to cross react with several other commercially important avian species and as such is a specific marker of B-cells for research in domestic galliforms. Clone 2C4 is unique in that it can be used in formalin fixed paraffin embedded tissues and is therefore of particular use on sections from infection studies.
Flow Cytometry	Use 10µl of the suggested working dilution to label 1x10 ⁶ cells in 100ul
Further Reading	1. Reddy, S.K. <i>et al.</i> (2008) The BAFF-Interacting receptors of chickens. Dev Comp Immunol. 32 (9): 1076-87. 2. Kothlow, S. <i>et al.</i> (2010) Prolonged effect of BAFF on chicken B cell development revealed by RCAS retroviral gene transfer in vivo. Mol Immunol. 47 (7-8): 1619-28. 3. Schneider, K. <i>et al.</i> (2004) Chicken BAFF--a highly conserved cytokine that mediates B cell survival. Int Immunol. 16 (1): 139-48. 4. Koskela, K. <i>et al.</i> (2004) Chicken B-cell-activating factor: regulator of B-cell survival in the bursa of fabricius. Scand J Immunol. 59 (5): 449-57.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	18 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10040 available at: 10040: https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP

Human Anti Mouse IgG1 (HCA036...) [HRP](#)

Recommended Useful Reagents

[MOUSE ANTI CHICKEN Bu-1a:RPE \(MCA2170PE\)](#)

[MOUSE ANTI CHICKEN Bu-1a:FITC \(MCA2170F\)](#)

[MOUSE ANTI CHICKEN Bu-1a/b \(MCA5764\)](#)

[MOUSE ANTI GUINEA FOWL Bu-1a/b \(MCA5958GA\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M286959:160520'

Printed on 01 May 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)