

Datasheet: MCA6020

Description:	MOUSE ANTI mCHERRY
Specificity:	mCHERRY
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	Clone 1
Isotype:	IgG1
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Western Blotting	■			10 ug/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	50% Glycerol
Approx. Protein Concentrations	IgG concentration 0.5 mg/ml
Immunogen	Recombinant mCherry expressed in <i>E. coli</i>
Specificity	Mouse anti mCherry antibody recognizes the red fluorescent mCherry protein, a monomeric engineered derivative of a fluorescent protein isolated from members of the Discosoma (mushroom coral) family. mCherry is a second generation fluorescent protein and by far the most popular red emitting fluorescent protein due to its high photostability and size (the other photostable red fluorescent protein tdTomato due to its tandem dimer nature is significantly larger than the 29 kDa mCherry protein). mCherry is optimally excited at a wavelength of 587nm and emits light at 610nm. Fluorescent proteins are used as tags in the research environment and are fused to the N- or C-terminus of proteins of interest. They are tools enabling the visualization of proteins, protein-protein interactions and are used in cell-sorting and FRET, FRAP and FLIM imaging experiments.

Western Blotting	Mouse anti mCherry antibody recognises a band of approximately 27 kDa on recombinant mCherry
Storage	Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10106 available at: 10106: https://www.bio-rad-antibodies.com/uploads/MSDS/10106.pdf
Regulatory	For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG1 (HCA036...)	HRP

Recommended Useful Reagents

[RABBIT ANTI GREEN FLUORESCENT PROTEIN \(AHP2984\)](#)
[RABBIT ANTI BLUE FLUORESCENT PROTEIN \(AHP2985\)](#)
[RABBIT ANTI CYAN FLUORESCENT PROTEIN \(AHP2986\)](#)
[RABBIT ANTI RED FLUORESCENT PROTEIN \(AHP2987\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M293660:160901'

Printed on 01 May 2018