

Datasheet: MCA6008GA

Description:	MOUSE ANTI SHEEP CD230 (aa143-154)
Specificity:	CD230 (aa143-154)
Other names:	PRP
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	ROS-BH1
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin	▪			
ELISA			▪	
Immunoprecipitation			▪	
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Sheep
Species Cross Reactivity	Reacts with: Bovine, Human, Hamster, Mouse, Red deer N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	Ig concentration 1.0 mg/ml
Immunogen	A truncated form of recombinant Ovine PRP spanning residues aa94 - aa233.

**External Database
Links**

UniProt:

[P23907](#) [Related reagents](#)
[P10279](#) [Related reagents](#)

Entrez Gene:

[493887](#) PRNP [Related reagents](#)
[281427](#) PRNP [Related reagents](#)

Synonyms

PRP, SIP

Fusion Partners

Spleen cells from immunised 129/Ola mice were fused with cells of the SP2/0 myeloma cell line.

Specificity

Mouse anti Sheep CD230 (Prpsc) antibody, clone ROS-BH1 recognizes the prion disease form of prion protein known as CD230. Transmissible spongiform encephalopathies (TSEs) or prion diseases are fatal infectious neurodegenerative diseases of humans and animals, and known as Scrapie (sheep and goats), BSE (cattle), CWD (cervine/deer) and CJD (humans). These diseases are biologically unique, as they are believed by some to be transmitted by an infectious agent comprised only of protein, with no nucleic acid component. Clinically, these diseases present with motor disturbances and behavioral changes. The major pathological changes seen are neuronal loss, vacuolation (spongiform change), proliferation and branching of glial cells, astrocytic proliferation and accumulation of the prion protein PrP^{Sc}, which can form amyloid plaques.

CD230. Also known as the prion protein (PrP) exists in two alternate forms; a normal cellular form (PrP^c) and a disease-associated form (PrP^{Sc}). The normal and pathological forms of the prion protein have identical amino acid sequences and differ only in their folded tertiary structure and biochemical properties.

Storage

Store at +4°C or at -20°C if preferred.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted.

Shelf Life

18 months from date of despatch.

Acknowledgements

This reagent was licenced from The Roslin Institute at the University of Edinburgh as part of the immunological toolbox project

**Health And Safety
Information**

Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

Rabbit Anti Mouse IgG (STAR8...) [DyLight®800](#)
Goat Anti Mouse IgG (STAR70...) [FITC](#)
Rabbit Anti Mouse IgG (STAR13...) [HRP](#)
Human Anti Mouse IgG1 (HCA036...) [HRP](#)
Goat Anti Mouse IgG (H/L) (STAR117...) [Alk. Phos.](#), [DyLight®488](#), [DyLight®549](#),
[DyLight®649](#), [DyLight®680](#), [DyLight®800](#),
[FITC](#), [HRP](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

'M318009:180718'

Printed on 01 Aug 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)