

Datasheet: MCA896GA

Description:	MOUSE ANTI SHEEP CD45
Specificity:	CD45
Other names:	LCA
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	VPM18
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/25 - 1/200
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA			■	
Immunoprecipitation	■			
Western Blotting			■	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Sheep
Species Cross Reactivity	Reacts with: Bovine Does not react with: Human N.B. Antibody reactivity and working conditions may vary between species.
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide (NaN ₃)
Carrier Free	Yes
Approx. Protein	IgG concentration 1.0 mg/ml

Concentrations

Immunogen Sheep T lymphocytes.

Fusion Partners Spleen cells from immunised BALB/c mice were fused with cells of the mouse NS0 myeloma cell line.

Specificity **Mouse anti Sheep CD45, clone VPM18** recognizes the ovine homologue of human CD45 otherwise known as LCA (leucocyte common antigen), a type I transmembrane glycoprotein and member of the protein-tyrosine phosphatase (PTP) family, expressed at high levels on all nucleated cells of haematopoietic origin, and which exists in several isoforms generated by the alternative splicing.

CD45 is a signal transduction protein vital for T and B cell antigen receptor-mediated activation, and is involved in many aspects of cellular processes, including growth, differentiation and oncogenic transformation. Studies have reported associations between CD45 and many cell surface antigens, including CD2, CD3, CD4, CD26 and CD90.

Clone VPM18 recognizes all isoforms of the CD45 molecule.

Flow Cytometry Use 10ul of the suggested working dilution to label 10^6 cells or 100ul whole blood

References

1. Hopkins, J. & Dutia, B.M. (1990) Monoclonal antibodies to the sheep analogues of human CD45 (leucocyte common antigen), MHC class I and CD5. Differential expression after lymphocyte activation in vivo. [Vet Immunol Immunopathol. 24 \(4\): 331-46.](#)
2. Pastoret P-P., *et. al.* (1998) Immunology of cattle (in Handbook of Veterbrate Immunology) Chapter XIII, p445

Storage Store at +4°C for one month or at -20°C for longer.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 18 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...)	RPE
Goat Anti Mouse IgG IgA IgM (STAR87...)	Alk. Phos. , HRP
Rabbit Anti Mouse IgG (STAR9...)	FITC
Goat Anti Mouse IgG (STAR77...)	HRP
Rabbit Anti Mouse IgG (STAR12...)	RPE
Goat Anti Mouse IgG (Fc) (STAR120...)	FITC , HRP
Rabbit Anti Mouse IgG (STAR8...)	DyLight®800

Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG1 (HCA036...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Useful Reagents

[MOUSE ANTI SHEEP CD45R:RPE \(MCA2221PE\)](#)

[MOUSE ANTI BOVINE CD14:FITC \(MCA2678F\)](#)

[MOUSE ANTI BOVINE CD14 \(MCA2678GA\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M319236:180719'

Printed on 27 Jul 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)