

Datasheet: MCA938GA

Description:	MOUSE ANTI RAT CD8 BETA
Specificity:	CD8 BETA
Format:	Purified
Product Type:	Monoclonal Antibody
Clone:	341
Isotype:	IgG1
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			1/100 - 1/200
Immunohistology - Frozen			■	
Immunohistology - Paraffin			■	
ELISA			■	
Immunoprecipitation	■			
Immunofluorescence	■			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	Purified IgG - liquid
Preparation	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Carrier Free	Yes
Approx. Protein Concentrations	IgG concentration 1 mg/ml
Immunogen	Rat/mouse T cell hybrids expressing CD8.
External Database	UniProt:

Links

[P05541](#) [Related reagents](#)

Entrez Gene:

[24931](#) Cd8b [Related reagents](#)

Synonyms

Cd8b1

Fusion Partners

Spleen cells from immunized BALB/c mice were fused with cells of the mouse X63Ag8.653 myeloma cell line.

Specificity

Mouse anti Rat CD8 beta antibody, clone 341 recognizes the beta chain of the rat CD8 cell surface antigen, expressed on thymocytes and peripheral cytotoxic/suppressor cells. CD8b is not expressed by NK cells or by intestinal intraepithelial T cells. Mouse anti Rat CD8 beta antibody, clone 341 is reported to block proliferative and cytotoxic responses of CD8 +ve effector cells *in vitro*.

Flow Cytometry

Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

References

1. Pérez-Cano, F.J. *et al.* (2010) Influence of breast milk polyamines on suckling rat immune system maturation. [Dev Comp Immunol. 34: 210-8.](#)
2. Torres-nagel, N. *et al.* (1992) Differential thymus dependence of rat CD8 isoform expression. [Eur J Immunol. 22 \(11\): 2841-8.](#)
3. Yan, Y. *et al.* (2003) Pathogenesis of Autoimmunity After Xenogeneic Thymus Transplantation [J Immunol. 170: 5936-46.](#)
4. Jörens, A. *et al.* (2010) Diabetes prevention by immunomodulatory FTY720 treatment in the LEW.1AR1-iddm rat despite immune cell activation. [Endocrinology. 151: 3555-65.](#)
5. Feferman, T. *et al.* (2005) Overexpression of IFN-induced protein 10 and its receptor CXCR3 in myasthenia gravis. [J Immunol. 174 \(9\): 5324-31.](#)

Storage

Store at +4°C for one month or at -20°C for longer.

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life

18 months from date of despatch.

Health And Safety Information

Material Safety Datasheet documentation #10040 available at: 10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory

For research purposes only

Related Products

Recommended Secondary Antibodies

Goat Anti Mouse IgG (STAR76...) [RPE](#)
Goat Anti Mouse IgG IgA IgM (STAR87...) [Alk. Phos.](#), [HRP](#)
Rabbit Anti Mouse IgG (STAR9...) [FITC](#)
Goat Anti Mouse IgG (STAR77...) [HRP](#)
Rabbit Anti Mouse IgG (STAR12...) [RPE](#)
Goat Anti Mouse IgG (Fc) (STAR120...) [FITC](#), [HRP](#)

Rabbit Anti Mouse IgG (STAR8...)	DyLight®800
Goat Anti Mouse IgG (STAR70...)	FITC
Rabbit Anti Mouse IgG (STAR13...)	HRP
Human Anti Mouse IgG1 (HCA036...)	HRP
Goat Anti Mouse IgG (H/L) (STAR117...)	Alk. Phos. , DyLight®488 , DyLight®549 , DyLight®649 , DyLight®680 , DyLight®800 , FITC , HRP

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL \(MCA1209\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M319237:180719'

Printed on 01 Aug 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)