

Datasheet: MCA938PE

Description:	MOUSE ANTI RAT CD8 BETA:RPE
Specificity:	CD8 BETA
Format:	RPE
Product Type:	Monoclonal Antibody
Clone:	341
Isotype:	IgG1
Quantity:	100 TESTS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	■			Neat

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Rat		
Product Form	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized		
Reconstitution	Reconstitute with 1 ml distilled water		
Max Ex/Em	Fluorophore	Excitation Max (nm)	Emission Max (nm)
	RPE 488nm laser	496	578
Preparation	Purified IgG prepared by affinity chromatography on Protein G from		
Buffer Solution	Phosphate buffered saline		
Preservative	0.09% Sodium Azide		
Stabilisers	1%	Bovine Serum Albumin	
	5%	Sucrose	
Immunogen	Rat/mouse T cell hybrids expressing CD8.		
External Database Links	UniProt: P05541 Related reagents Entrez Gene:		

Synonyms	Cd8b1
Fusion Partners	Spleen cells from immunized BALB/c mice were fused with cells of the mouse X63Ag8.653 myeloma cell line.
Specificity	Mouse anti Rat CD8 beta antibody, clone 341 recognizes the beta chain of the rat CD8 cell surface antigen, expressed on thymocytes and peripheral cytotoxic/suppressor cells. CD8b is not expressed by NK cells or by intestinal intraepithelial T cells. Mouse anti Rat CD8 beta antibody, clone 341 is reported to block proliferative and cytotoxic responses of CD8 +ve effector cells <i>in vitro</i> .
Flow Cytometry	Use 10ul of the suggested working dilution to label 10 ⁶ cells in 100ul.
References	1. Pérez-Cano, F.J. <i>et al.</i> (2010) Influence of breast milk polyamines on suckling rat immune system maturation. Dev Comp Immunol. 34: 210-8. 2. Torres-nagel, N. <i>et al.</i> (1992) Differential thymus dependence of rat CD8 isoform expression. Eur J Immunol. 22 (11): 2841-8. 3. Yan, Y. <i>et al.</i> (2003) Pathogenesis of Autoimmunity After Xenogeneic Thymus Transplantation J Immunol. 170: 5936-46. 4. Jörns, A. <i>et al.</i> (2010) Diabetes prevention by immunomodulatory FTY720 treatment in the LEW.1AR1-iddm rat despite immune cell activation. Endocrinology. 151: 3555-65. 5. Feferman, T. <i>et al.</i> (2005) Overexpression of IFN-induced protein 10 and its receptor CXCR3 in myasthenia gravis. J Immunol. 174 (9): 5324-31.
Storage	Prior to reconstitution store at +4°C. Following reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet documentation #10075 available at: 10075: https://www.bio-rad-antibodies.com/uploads/MSDS/10075.pdf
Regulatory	For research purposes only

Related Products

Recommended Negative Controls

[MOUSE IgG1 NEGATIVE CONTROL:RPE \(MCA1209PE\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M305239:170426'

Printed on 01 Aug 2018

