

Datasheet: OBT1530

Description:	MOUSE ANTI STAPHYLOCOCCUS AUREUS
Specificity:	STAPHYLOCOCCUS AUREUS
Format:	Ascites
Product Type:	Monoclonal Antibody
Clone:	11-248.2
Isotype:	IgM
Quantity:	0.1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			▪	
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1/16,000
Immunoprecipitation			▪	
Western Blotting			▪	

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Bacterial
Product Form	Ascites fluid
Preservative Stabilisers	0.1% Sodium Azide (NaN ₃)
Immunogen	UV inactivated <i>Staphylococcus aureus</i> cells.
Specificity	Mouse anti <i>Staphylococcus aureus</i> antibody, clone STAPH 11-248.2 recognizes <i>Staphylococcus aureus</i>. No reactivity is seen with <i>Streptococcus</i> sp., <i>Mycoplasma</i> sp., or <i>Corynebacteria</i> sp.
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 12 months from date of despatch.

Health And Safety Information Material Safety Datasheet Documentation #10342 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10342.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M285886:160414'

Printed on 25 May 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)