

Datasheet: OBT1508P

Description:	RABBIT ANTI MOUSE IgG1:HRP
Specificity:	IgG1
Format:	HRP
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry			■	
Immunohistology - Frozen	■			
Immunohistology - Paraffin			■	
ELISA	■			1/20000 - 1/100000
Western Blotting	■			1/2000 - 1/10000
Immunoblotting	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilised
Reconstitution	Reconstitute with 1.0 ml distilled water
Preparation	Purified IgG prepared by affinity chromatography
Buffer Solution	Phosphate buffered saline
Preservative	0.01% gentamicin sulphate
Stabilisers	1% Bovine serum albumin
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
Immunogen	Mouse IgG1 gamma 1 heavy chain

External Database Links

UniProt:

[P01869](https://www.uniprot.org/entry/P01869)

[Related reagents](#)

Entrez Gene:

[16017](#) IgHg1 [Related reagents](#)

[16017](#) IgHg1 [Related reagents](#)

Synonyms	Igh-4
Specificity	Rabbit anti Mouse IgG1 antibody is specific for mouse IgG1 and shows minimal (<0.5%) cross-reactivity with mouse IgG2a, IgG2b, IgG3 and IgM.
References	1. Kerr K <i>et al.</i> (2010) Inflammatory cytokine responses in a pregnant mouse model of <i>Chlamydomphila abortus</i> infection. Vet Microbiol. 144 (3-4): 392-8.
Storage	Prior to reconstitution store at +4°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	12 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet documentation available at: Material Safety Datasheet Documentation #10482 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10482.pdf
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[AbGUARD® HRP STABILIZER PLUS \(BUF052A\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052B\)](#)

[AbGUARD® HRP STABILIZER PLUS \(BUF052C\)](#)

[TMB CORE \(BUF056A\)](#)

[TMB CORE+ \(BUF062A\)](#)

[TMB SIGNAL+ \(BUF054A\)](#)

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M314507:180412'

Printed on 13 Apr 2018