

Datasheet: PEP003

Description:	RECOMBINANT HORSE TNF ALPHA
Name:	TNF ALPHA
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	20 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA			■	
Functional Assays	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Horse
Product Form	Purified Recombinant Protein- lyophilised
Reconstitution	Reconstitute with 0.5 ml distilled water. Further dilutions should be made in a buffer containing carrier protein. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Note	Further dilutions should be made in a buffer containing carrier protein. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>pichia pastoris</i>
Preservative Stabilisers	None present
Endotoxin Level	0.15EU/ug
External Database Links	UniProt: P29553 Related reagents

Entrez Gene:[100033834](#) TNF [Related reagents](#)

Synonyms TNFA, TNFSF2

Product Information **Equine Tumour Necrosis Factor Alpha** (TNF α) is produced by a range of cells, including macrophages and monocytes and is an important mediator of the equine inflammatory response.

Equine Tumour Necrosis Factor Alpha has an ED₅₀ < 0.05ng/ml, determined by measuring the cytolytic effect of recombinant equine TNF α on murine WEHI-13VAR-164 cells in the presence of actinomycin D.

Protein Molecular Weight 17.1kDa

Purity >95% by SDS PAGE analysis

Storage Prior to reconstitution store at +4°C.
Following reconstitution store at -70°C.
Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 12 months from date of reconstitution.

Health And Safety Information Material Safety Datasheet Documentation #10268 available at:
<https://www.bio-rad-antibodies.com/uploads/MSDS/10268.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M260188:140912'

Printed on 07 Jun 2017

© 2017 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)