

Datasheet: PHP045A

Description:	RECOMBINANT HUMAN INTERLEUKIN-6
Name:	IL-6
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.2 - 0.4ng/well
Western Blotting	▪			1.5 - 3.0ng/lane
Functional Assays	▪			0.1 - 10ng/ml

Where this protein has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the protein for use in their own system using appropriate positive/negative controls.

Target Species	Human
Product Form	Purified recombinant protein - lyophilised
Reconstitution	Reconstitute with 50ul distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For extended storage, the addition of 0.1% bovine serum albumin (BSA) is recommended.
Preparation	Recombinant protein expressed in <i>E. coli</i>
Preservative Stabilisers	None present
Endotoxin Level	<0.1 ng/ug
Approx. Protein Concentrations	0.1 mg/ml after reconstitution.
External Database Links	UniProt: P05231 Related reagents Entrez Gene: 3569 IL6 Related reagents

Synonyms	IFNB2
Product Information	Recombinant Human Interleukin-6 protein produced in <i>E. coli</i> constitutes full length non-glycosylated human IL-6 and retain full biological activity. Interleukin-6 is a potent lymphoid cell growth factor, acting on B cells, T cells and hybridoma cells. This preparation has an ED50 of 0.1ng/ml when assayed by proliferation of murine 7TD1 cells.
Protein Molecular Weight	20.9 kDa (184 amino acids)
Activity	>1 x 10 ⁷ units/mg
Purity	>98% by SDS PAGE analysis and HPLC analysis
ELISA	PHP045A may be used as the standard in ELISA applications with either AHP1040 or AHP1040B .
Western Blotting	PHP045A may be used as a positive control for western blotting applications with AHP1040 or AHP1040B .
References	1. Ban, S.A. <i>et al.</i> (2014) Combined immunodeficiency evolving into predominant CD4+ lymphopenia caused by somatic chimerism in JAK3. J Clin Immunol. 34 (8): 941-53.
Storage	Prior to reconstitution store at +4°C. Following reconstitution store at -20°C. This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	3 months from date of reconstitution
Health And Safety Information	Material Safety Datasheet documentation #10527 available at: 10527: https://www.bio-rad-antibodies.com/uploads/MSDS/10527.pdf
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[RABBIT ANTI HUMAN INTERLEUKIN-6 \(AHP1040\)](#)
[RABBIT ANTI HUMAN INTERLEUKIN-6:Biotin \(AHP1040B\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M318210:180718'

Printed on 01 Aug 2018