

Datasheet: PMOP001

Description:	MONKEY IgG
Name:	IgG
Format:	Purified
Product Type:	Purified Protein
Isotype:	IgG
Quantity:	2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Monkey
Product Form	Purified IgG - lyophilised
Reconstitution	Reconstitute with 1.0 ml distilled water immediately prior to use. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified monkey IgG was prepared from normal serum by ion exchange chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.01% Sodium Azide (NaN ₃)
Approx. Protein Concentrations	IgG concentration 2.0mg/ml
Product Information	PMOP001 is a preparation of purified polyclonal monkey IgG prepared from normal serum. By immunoelectrophoresis, PMOP001 was shown to react with anti monkey IgG and anti-monkey serum.
Storage	Prior to reconstitution store at +4°C. After reconstitution store at -20°C.

Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life 18 months from date of reconstitution

Health And Safety Information Material Safety Datasheet documentation #10258 available at:
10258: <https://www.bio-rad-antibodies.com/uploads/MSDS/10258.pdf>

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M317244:180619'

Printed on 21 Jun 2018

© 2018 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)