

Datasheet: PMP38

Description:	RECOMBINANT MOUSE INTERLEUKIN-2
Name:	IL-2
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	20 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.2 - 0.4 ng/well
Western Blotting	▪			1.5 - 3.0 ng/lane
Functional Assays	▪			0.1 - 1.0 ng/ml

Where this protein has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate positive/negative controls.

Target Species	Mouse
Product Form	Purified recombinant protein - lyophilised
Reconstitution	Reconstitute with 0.2ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>E.coli</i>
Buffer Solution	10mM Sodium Citrate pH4.0
Preservative Stabilisers	None present
Endotoxin Level	<0.1 ng/ug
Approx. Protein Concentrations	0.1 mg/ml after reconstitution.
External Database Links	UniProt: P04351 Related reagents Entrez Gene:

Synonyms	IL-2
Product Information	Interleukin-2 (IL-2) is a potent lymphoid cell growth factor which primarily acts on T-cells.
Protein Molecular Weight	17.2kD (149 amino acid sequence)
Activity	5 x 10 ⁶ Units/mg
Purity	>98% by SDS PAGE and HPLC
ELISA	This product may be used as a standard in ELISA applications either AAM35 or AAM35B .
Western Blotting	This product may be used as the positive control in Western Blotting applications either AAM35 or AAM35B .
References	1. Bénéteau, M. <i>et al.</i> (2012) Combination of glycolysis inhibition with chemotherapy results in an antitumor immune response. Proc Natl Acad Sci U S A. 109: 20071-6.
Storage	Prior to reconstitution store at -20°C. Following reconstitution store at -20°C. This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	3 months from date of reconstitution. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10175 available at: 10175: https://www.bio-rad-antibodies.com/uploads/MSDS/10175.pdf
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M316904:180614'

Printed on 21 Jun 2018