

Datasheet: PMP47

Description:	RECOMBINANT MOUSE INTERLEUKIN-15
Name:	IL-15
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA			■	
Functional Assays	■			0.1ng/ml - 1.0ng/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified recombinant protein - lyophilised
Reconstitution	Reconstitute with 10 ul distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For extended storage the addition of 0.1% BSA is recommended.
Preparation	Recombinant protein expressed in <i>E. coli</i>
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Endotoxin Level	<0.1 ng/ug
Approx. Protein Concentrations	1.0 mg/ml
External Database Links	UniProt: P48346 Related reagents Entrez Gene:

Product Information	Mouse Interleukin 15 (IL-15) is a novel cytokine which affects T-cell activation and proliferation.
----------------------------	--

Protein Molecular Weight	13.3kDa (115 amino acid sequence)
---------------------------------	-----------------------------------

Activity	2 x 10 ⁵ Units/mg
-----------------	------------------------------

Purity	>98% by SDS PAGE and HPLC
---------------	---------------------------

Storage	Prior to reconstitution store at +4°C. Following reconstitution store at -20°C.
----------------	---

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life	3 months from date of reconstitution.
-------------------	---------------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10527 available at: 10527: https://www.bio-rad-antibodies.com/uploads/MSDS/10527.pdf
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

'M318290:180718'

Printed on 01 Aug 2018