

Datasheet: PMP77

Description:	RECOMBINANT MOUSE SDF-1 ALPHA
Name:	SDF-1 ALPHA
Other names:	CXCL12a
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified recombinant protein - lyophilised
Reconstitution	Reconstitute with 100ul distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified recombinant SDF-1 alpha expressed in <i>E.coli</i>
Preservative Stabilisers	None present
Endotoxin Level	<1EU/ug
Approx. Protein Concentrations	0.1mg/ml
External Database Links	UniProt: P40224 Related reagents Entrez Gene: 20315 Cxcl12 Related reagents
Synonyms	Sdf1

Product Information	SDF-1 alpha, otherwise known as CXCL12a, is a stromal cell derived CXC chemokine which, together with SDF-1 beta, arises from alternative splicing of the SDF-1 gene, and acts as a chemoattractant for T and B lymphocytes, monocytes and migratory neurons. SDF-1 binds with high-affinity to the G protein-coupled receptor CXCR4 (fusin), and acts as an inhibitor of the CXCR4-mediated entry of HIV-1 virus into target T-cells.
Molecular Weight	7.9kD (68 amino acid residue homodimer)
Activity	Determined by its ability to chemoattract human peripheral blood monocytes, using a concentration of 50-100 ng/ml.
Purity	>98% by SDS PAGE and HPLC analysis
Storage	Prior to reconstitution store at +4°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Shelf Life	3 months from date of reconstitution.
Health And Safety Information	Material Safety Datasheet Documentation #10268 available at: https://www.bio-rad-antibodies.com/uploads/MSDS/10268.pdf
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

'M229878:130108'

Printed on 25 May 2017